

CÁPSULAS Y NAVECILLAS PARA ANALIZADORES ELEMENTALES

[click
aqui](#)

PRODUCTOS PARA ANALIZADORES ELEMENTALES LECO

[click
aqui](#)

PRODUCTOS PARA ANALIZADORES ELEMENTALES ELEMENTAR

[click
aqui](#)

PRODUCTOS PARA ANALIZADORES ELEMENTALES THERMO

[click
aqui](#)

PRODUCTOS PARA ANALIZADORES ELEMENTALES EUROVECTOR

[click
aqui](#)

MATERIALES DE REFERENCIA CERTIFICADOS / PATRONES

[click
aqui](#)

PATRONES / STANDARDS LECO

[click
aqui](#)

COBRE PARA ANALIZADORES ELEMENTALES

[click
aqui](#)

PRODUCTOS PARA ICP PERKIN ELMER

[click
aqui](#)





CÁPSULAS Y NAVECILLAS PARA ANALIZADORES ELEMENTALES

Listado de Artículos



EQUILAB

SANTIS
analytical

Todas nuestras cápsulas están especialmente fabricadas a semejanza de las originales y son rigurosamente controladas para cumplir los estándares de calidad, garantizar bajos valores de blancos, alta calidad y peso uniforme. Debido a la buena cooperación con los principales fabricantes de instrumentos analíticos nuestros productos son siempre de tecnología de última generación.

Podemos fabricar cápsulas, discos y navecillas en tamaños y materiales diferentes a los que aparecen en la lista. Así que por favor no dude en contactar con nosotros.



Article List for Pressed Capsules, Smooth Wall Capsules and Weighing Boats

All our capsules are manufactured to exacting standards and are quality controlled to ensure consistent low blank values, high quality and uniform weight. Due to the good co-operation with the leading manufacturers of analytical instruments our products are always of the latest state-of-the-art.

We can manufacture capsules, discs and boats in sizes and materials other than those listed above. So please don't hesitate to contact us.

Pressed Capsules for solid samples

Precleaned capsules are carefully solvent extracted then individually packaged in trays of 96 capsules ready for immediate use.

Säntis No.	Description	Diameter mm	Height mm	Volume ml	Qty./ pack.
Pressed Tin Capsules					
SA76980002	Tin capsules for solids	3.20	4.00	0.020	250
SA76980202	Tin capsules for solids	3.20	6.00	0.040	250
SA76980501	Tin capsules for solids	3.30	5.00	0.030	100
SA76980502	Tin capsules for solids	3.30	5.00	0.030	250
SA76980702	Tin capsules for solids	4.00	6.00	0.070	250
SA76980703	Tin capsules for solids / Precleaned	4.00	6.00	0.070	96
SA76980901	Tin capsules for solids	5.00	8.00	0.120	100
SA76980902	Tin capsules for solids	5.00	8.00	0.120	250
SA76980903	Tin capsules for solids / Precleaned	5.00	8.00	0.120	96
SA76981101	Tin capsules for solids	5.00	9.00	0.150	100
SA76981102	Tin capsules for solids	5.00	9.00	0.150	250
SA76981103	Tin capsules for solids	5.00	9.00	0.150	400
SA76981302	Tin capsules for solids	5.00	12.00	0.220	250
SA76981501	Tin capsules for solids	8.00	11.00	0.550	100
SA76156701	Tin capsules for solids	6.70	14.00	0.500	50
SA999926	Tin capsules for solids	8.00	20.00	0.950	100
SA76162401	Tin capsules for solids	10.00	10.00	0.800	100
SA76152301	Tin foil Cups				100
SA502397400	Large Tin Foil Cups				400

Pressed Tin Capsules - Light (weigh 33% less than the standard)

SA76980002L	Tin capsules for solids / light	3.20	4.00	0.020	250
SA76980502L	Tin capsules for solids / light	3.30	5.00	0.030	250
SA76980702L	Tin capsules for solids / light	4.00	6.00	0.070	250
SA76980703L	Tin capsules for solids / light / Precleaned	4.00	6.00	0.070	96
SA76980901L	Tin capsules for solids / light	5.00	8.00	0.120	100
SA76980902L	Tin capsules for solids, light	5.00	8.00	0.120	250
SA76980903L	Tin capsules for solids / light / Precleaned	5.00	8.00	0.120	96
SA76981101L	Tin capsules for solids / light	5.00	9.00	0.150	100
SA76981102L	Tin capsules for solids / light	5.00	9.00	0.150	250
SA76981401L	Tin capsules for solids / light	4.00	4.75	0.050	100
SA76981601L	Tin capsules for solids / light	7.00	11.50	0.350	100

Tin Discs

Säntis No.	Description	Diameter mm	Height mm	Volume ml	Qty./ pack.
SA999917	Tin discs, thickness 0.010 mm	21.00	-	-	100
SA999917L	Tin discs / light, thickness 0.008 mm	21.00	-	-	100
SA999901	Tin discs, thickness 0.010 mm	28.00	-	-	100
SA999918	Tin discs, thickness 0.010 mm	30.00	-	-	100
SA999918L	Tin discs / light, thickness 0.008 mm	30.00	-	-	100
SA999902	Tin discs, thickness 0.010 mm	35.00	-	-	100
SA999972	Tin discs, thickness 0.010 mm	53.00	-	-	100
SA999976	Tin discs, thickness 0.010 mm	62.00	-	-	100
SA999976L	Tin discs / light, thickness 0.008 mm	65.00	-	-	100
SA990803	Tin discs, thickness 0.010 mm	86.00	-	-	100

Tin Squares

SA999990	Tin foil Squares 35 x 35 mm	35.00	-	-	400
SA999991	Tin foil Squares 50 x 50 mm	50.00	-	-	450

Pressed Silver Capsules

SA76980006	Silver capsules for solids	3.20	4.00	0.020	250
SA76980206	Silver capsules for solids	3.20	6.00	0.040	250
SA76980505	Silver capsules for solids	3.30	5.00	0.030	100
SA76980506	Silver capsules for solids	3.30	5.00	0.030	250
SA76980705	Silver capsules for solids	4.00	6.00	0.070	100
SA76980706	Silver capsules for solids	4.00	6.00	0.070	250
SA76980906	Silver capsules for solids	5.00	8.00	0.120	250
SA76981105	Silver capsules for solids	5.00	9.00	0.150	100
SA76981106	Silver capsules for solids	5.00	9.00	0.150	250
SA76981306	Silver capsules for solids	5.00	12.00	0.220	250
SA76162406	Silver capsules for solids	10.00	10.00	0.800	100

Silver Squares

SA999992	Silver foil Squares 40 x 40 mm	40.00	-	-	100
----------	--------------------------------	-------	---	---	-----

Smooth wall capsules for liquid samples

Flat Bottom Smooth Tin Capsules

SA84990000	Tin capsules for liquids	1.06	5.50	0.003	50
SA84990002	Tin capsules for liquids	1.06	5.50	0.003	250
SA84990300	Tin capsules for liquids	2.00	5.00	0.012	50
SA84990301	Tin capsules for liquids	2.00	5.00	0.012	100
SA84990302	Tin capsules for liquids	2.00	5.00	0.012	250
SA84990600	Tin capsules for liquids	2.90	6.00	0.025	50
SA84990601	Tin capsules for liquids	2.90	6.00	0.025	100
SA84990602	Tin capsules for liquids	2.90	6.00	0.025	250
SA84991201	Tin capsules for liquids	3.50	4.00	0.030	100
SA84991202	Tin capsules for liquids	3.50	4.00	0.030	250
SA84991502	Tin capsules for liquids	3.50	5.50	0.040	250

Säntis No.	Description	Diameter mm	Height mm	Volume ml	Qty./ pack.
SA84991802	Tin capsules for liquids	3.50	7.00	0.060	250
SA84992100	Tin capsules for liquids	3.50	9.00	0.070	50
SA84992101	Tin capsules for liquids	3.50	9.00	0.070	100
SA84992102	Tin capsules for liquids	3.50	9.00	0.070	250
SA84992702	Tin capsules for liquids	4.75	11.00	0.150	250
SA84992801B	Tin capsules for liquids	5.00	12.00	0.180	100
SA84992801	Tin capsules for liquids	5.00	13.00	0.180	100
SA84992801C	Tin capsules for liquids	5.00	18.00	0.300	100
SA84992901	Tin capsules for liquids	6.00	16.00	0.300	100
SA84993001	Tin capsules for liquids	6.75	14.00	0.450	100
SA84994401	Tin capsules for liquids	8.00	20.00	0.980	100

Flat Bottom Smooth Silver Capsules

SA84990305	Silver capsules for liquids	2.00	5.00	0.021	50
SA84990306	Silver capsules for liquids	2.00	5.00	0.012	250
SA84991206	Silver capsules for liquids	3.50	4.00	0.025	250
SA84991506	Silver capsules for liquids	3.50	5.50	0.040	250
SA84992106	Silver capsules for liquids	3.50	9.00	0.070	250

Flat Bottom Smooth Nickel Capsules

SA76980801	Nickel capsules for liquids	4.50	10.00	0.050	100
SA76981701	Nickel capsules for liquids	7.00	16.00	0.046	100

Pressed Boats

Säntis No.	Description	Width mm	Length mm	Height mm	Qty./ pack.
------------	-------------	----------	-----------	-----------	-------------

Pressed Tin Boats

SA22137418	Tin Boats	4.00	4.00	11.00	1000
SA22137419	Tin Boats	6.00	6.00	12.00	500
SA22137420	Tin Boats	8.00	8.00	15.00	250
SA76982502	Tin Boats	6.50	3.50	2.50	250
SA76982802	Tin Boats	10.00	4.00	4.50	250
SA76983102	Tin Boats	12.00	4.00	2.70	250

Pressed Aluminium Boats

SA76982801	Aluminium Boats	10.00	4.00	4.50	250
SA76983402	Aluminium Boats	12.00	4.00	4.50	250
SA76983701	Aluminium Boats	25.00	7.00	7.00	50

Pressed Silver Boats

SA76982506	Silver Boat	6.50	3.50	2.50	250
SA22133212	Silver Boats	6.00	6.00	12.00	100
SA22133213	Silver Boats	8.00	8.00	16.00	70
SA22133214	Silver Boats	10.00	10.00	20.00	42

Säntis No.	Description	Diameter mm	Height mm	Volume ml	Qty./ pack.
------------	-------------	----------------	--------------	--------------	----------------

Weighing Pans

Weighing Pans

SA76985002	Aluminium Weighing Pan	6.00	2.50		250
SA76985301	Aluminium Weighing Pan	13.00	3.50		100



Leco

Productos para Analizadores Elementales

Listado de Artículos



Article List for Elemental Analyzer

Leco™

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
145-601	SAAR5601	Water Temp Regulator	
156-002	SAAR002	"O" Ring (Upper Electrode)	
156-004	SAAR005	Teflon Backup Seal	
163-012	SAAR63012	Tubing Insert	
170-001	SAAR0001	Spring, Lower Electrode	
172-505	SAAR2505	Teflon Tubing, 1/4" OD	
175-574	SAAR5574	Optical Switch Assembly	
175-590	SAAR5590	Pump Repair Kit, SC32/132	
175-631	SAAR5631	Motor Brush for 774-571	
175-632	SAAR5632	"O" Ring, Flow Controller	
175-659	SAAR5659	Paper Deflector V55147	
175-674	SAAR5674	Tear Off Blade Spring	
188-012	SAAR293	Small Throat Brush	
190-142	SAAR142	Screws for Lower Electrode, 4 Required	
190-739	SAAR739	Nylon Spacer Screw	
194-060	SAAR4060	Hex Jam Nut for AR5814	
204-001	SAAR001	Thermoregulator TC30, TN14n R016	
501-001	SAAR1701	Coal Standard, 0.4% Sulfur	
501-005	SAAR1704	Coal Standard, 1.0% Sulfur	
501-008	SAAR008	Alphacel II Tungsten	2.2kg
501-009	SAAR1009	Spring Clip for Dust Trap	
501-010	SAAR010	Dust Trap Complete	
501-011	SAAR011	Dust Trap Body Insert	
501-012	SAAR011B	Dust Trap Reservoir	
501-013	SAAR013	Filter Cloths	10 pcs
501-014	SAAR1705	Coal Standard, 1.5% Sulfur	
501-020	SAAR1708	Coal Standard, 3.0% Sulfur	
501-022	SAAR1709	Coal Standard, 3.5% Sulfur	
501-024	SAAR322	Cast Iron Standard 3.41% C, 0.050% S	250g
501-027	SAAR204	Aspirator Bulb	
501-031	SAAR031	Induction Coil, Hydrocarbon	
501-034	SAAR1034	Calcium Carbonate, 12%	50g
501-045	SAAR1045	Zircon Crucible, CHN 900	1 pc.
501-046	SAAR2021	Sucrose Standard	50g
501-053	SA990744	Acetanilide Standard	1g
501-057	SAAR079G	Alphanal	3.785kg
501-059	SA84992801	Tin Capsules, 5x13mm	100 pcs.
501-060	SAAR060	Manganese Dioxide 30 CC	
501-061	SAAR1061	Arrowroot Starch	227g
501-062	SAAR061	Boat Puller	
501-067	SAAR067	Rotometer Ball Only	
501-073	SAAR073	Graphite Powder (Bottle)	100g
501-076	SAAR076	Tin Metal	
501-077	SAAR077	Iron Chip	
501-078	SAAR078	Iron Powder	
501-079	SAAR079	Alphanal Wetting Agent	118g
501-081	SAAR081	Glass Wool	
501-082	SAAR022	Brush	
501-083	SAAR304	Ductile Iron Standard, 4.23% C	
501-089	SAAR507	Valve Assembly for 516-000	
501-095	SAAR2814	Residual Oil, 1.0% S	
501-099	SAAR2837	Lube Oil, 4% S	
501-101	SAAR800	Arrowroot Starch Solution	
501-105	SAAR301	Cast Iron Standard 2.02% C, 0.018% S	150g

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
501-110	SAAR110	3/8" I.D. Silicone Tubing	
501-110-019	SAAR11019	Silicone Rubber Tube, 3/8" x 19"	
501-112	SAAR112	1/4" I.D. Neoprene Tubing	
501-123	SAAR307	Tungsten Carbide Standard 6.15% C, 0.0041% S	100g
501-130	SAAR1710	Coal Standard, 4.0% Sulfur	
501-144	SAAR144	Magnesium Oxide	
501-147	SAAR148	1g Copper Pin Standard High Oxygen 489 PPM O, 6 PPM S	100 pcs.
501-148	SAAR149	1g Copper Pin Standard High Oxygen 323 PPM O, 11 PPM S	100 pcs.
501-165	SAAR1165	Drierite Indicating 10-20M	454g
501-170	SA501170	Rare Earth Copper Oxide	50g
NEW: aquatak® Replacement for Anhydrone			
	SA100050	aquatak®	50g
	SA100100	aquatak®	100g
	SA100500	aquatak®	500g
501-171	SAAR171	Anhydrone, 8-20 Mesh	454g
501-179	SAAR079P	Alphanal	454g
501-184	SAAR1711	Coal Standard, 5.0% Sulfur	
501-197	SAAR775	WR-12 Trap Material	
501-233	SAAR233	Ceramic Adhesive Alumina	57g
501-234	SAAR234	Oxygen Specimen File	
501-241	SAAR241	Silicone Grease	150g
501-263	SAAR263	Copper Metal Accelerator	1.36kg
501-264	SAAR264	Copper Metal Accelerator	227g
501-265	SAAR1265	Crucible Holder Tray	
501-272	SAAR2017	Corn Gluten SRM	
501-280	SAAR634	"O" Ring	
501-291	SAAR291	Oxygen Regulator Assembly	
501-369	SAAR1727	Prox Coal Standard 2.5% S	
501-426	SAAR426	Combustion Control, Solid Samples	
501-427	SAAR427	Combustion Control, Liquid Samples	
501-428	SAAR1721	Prox Coal Standard 0.4% S	
501-429	SAAR1723	Prox Coal Standard 1.0% S	
501-430	SAAR1724	Prox Coal Standard, 1.5% S	
501-431	SAAR1726	Prox Coal Standard, 2.0% S	
501-432	SAAR1728	Prox Coal Standard, 3.0% S	
501-433	SAAR1730	Prox Coal Standard 4.0% S	
501-438	SAAR734	Prox Coal Standard 0.9% S	
501-439	SAAR000	Parafin Oil, 118 mL/Btl C, H, N BTU	
501-441	SAAR2021	Sucrose	50g
501-457	SAAR400	Tin Coated Copper	1.5kg
501-458	SAAR2822	#2 Diesel Oil	
501-501	SAAR870	1g C/S Ring Standard 0.032% C, 0.033% S	454g
501-502	SAAR871	1g C/S Ring Standard 0.051% C, 0.034% S	454g
501-503	SAAR872	1g C/S Ring Standard 0.180% C, 0.019% S	454g
501-504	SAAR873	1g C/S Ring Standard 0.468% C, 0.0109% S	454g
501-505	SAAR874	1g C/S Ring Standard 0.723% C, 0.022% S	454g
501-506	SAAR875	1g C/S Ring Standard 0.829% C, 0.0067% S	454g
501-529	SAAR556	1g Steel Ball, 5.7ppm H2	
501-551	SAAR660	0.5g O/N in Steel Pin Standard 0.0045% O, 0.0060% N	100 pcs
501-552	SAAR656	0.5g O/N in Steel Pin Standard 0.0081% O, 0.0213% N	100 pcs
501-561	SAAR2016	Soy Bean Meal SRM	
501-563	SAAR2025	Corn Meal SRM	
501-571	SAAR571	Copper Capsules	100 pcs
501-587	SAAR912	Platinized Silica	15g

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
501-592	SAAR062	Manganese Dioxide, 250CC	
501-598	SAAR598	Nibbled Nickel Flux	100g
501-608	SAAR1608	Quartz Wool Wadding	200g
501-609	SAAR615	Furnace Regent	100g
501-614	SAAR1614	Spatula	
501-621	SA501621W	Copper Turnings	200g
501-632	SAAR632	Mole Sieve Trap	14g
501-636	SA990709	Vanadium Pentoxide	1g
501-640	SAAR366	Copper Metal Accelerator, 25 lb	
501-641	SAAR643	Copper Metal Scrubber Material	454g
501-644	SAAR644	1g O/N in Steel Pin Standard 3PPM O, 157PPM N	100 pcs.
501-645	SAAR1651	1g O/N in Steel Pin Standard 0.0057% O, 0.0210% N	100 pcs.
501-648	SAAR2776	Ultimate Coal Standard	50g
501-673	SAAR673	Wrought Iron Chip Standard 0.0004% C, 0.0017% S	454g
501-674	SAAR881	1g C/S Pin Standard 0.014% C, 0.0045% S	454g
501-675	SAAR882	1g C/S Pin Standard 0.063% C, 0.023% S	454g
501-676	SAAR890	1g C/S Pin Standard 0.099% C, 0.0182% S	454g
501-677	SAAR895	1g C/S Pin Standard 0.247% C, 0.021% S	454g
501-678	SAAR884	1g C/S Pin Standard 0.4410% C, 0.0110% S	454g
501-679	SAAR885	1g C/S Pin Standard 0.805% C, 0.022% S	454g
501-697	SAAR612	Platinized Charcoal	50g
501-951	SAAR1951	Nickel Calibration Wire	
501-955	SAAR1722	Prox Coal Standard .7% S	
501-957	SAAR1729	Prox Coal Standard 3.5% S	
501-961	SAAR304	Ductile Iron Standard 4.23% C, 0.014% S	150g
501-964	SAAR1712	Coal Standard 6.0% Sulfur	
501-967	SAAR2851	Extreme Pressure Oil	
501-968	SAAR2852	Extreme Pressure Oil	
501-969	SAAR2853	Extreme Pressure Oil	
502-005	SAAR1205	Pipette Brush	
502-007	SAAR2189	Tube Cleaner, 8mm x 4ft, FP428	
502-008	SA84993701	Tin Plugs	150 pcs.
502-016	SAAR655	0.5g O/N in Steel Pin Standard 0.0021% O, 0.393% N	100 pcs.
502-020	SAAR599	Etchable Nickel Flux	100g
502-023	SAAR473	Reagent Tube Funnel	
502-024	SAAR651	0.1g H/O/N in Titanium Standard 28 PPM H, 1320 PPM O, 49 PPM N	100 pcs.
502-026	SAAR2818	Residual Oil	
502-032	SAAR078	Accelerator, Fe Powder	454g
502-040	SA84992901	Tin Capsule, 6x16mm	100 pcs.
502-042	SAAR025	Quartz Crucible, CHN	
502-049	SAAR049	Pt-Calalyst	50g
502-060	SAAR556	1g Hydrogen in Steel Pin Standard 5.0 PPM H	100 pcs.
502-061	SAAR555	1g Hydrogen in Steel Pin Standard 3.4 PPM H	100 pcs.
502-062	SAAR4020	Soil Standard, 1% C	100g
502-071	SAAR2864	High Sulfur Residual Oil, 118mL 30% S	
502-072	SAAR655	0.5g O/N in Steel Pin Standard 0.0021% O, 0.393% N	100 pcs.
502-076	SAAR571	Copper Capsules	100 pcs.
502-077	SA84993701	Tin Plugs	150 pcs.
502-085	SAAR3085	Zn Sulfide, 32%	50g
502-091	SAAR091	Calcium Oxalate Standard	50g
502-092	SA502092	EDTA Standard	50g
502-122	SAAR2122	Firing Fuse Wire, 2"	10 pcs.
502-134	SAAR2134	Benzoic Acid Tablet 1.45g	10 pcs.
502-134	SAAR2134V	Benzoic Acid Tablet 1.45g	50 pcs.

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
502-154	SAAR648	0.1g H/O/N in Titanium Standard 0.0161% H, 0.158% O, 0.015% N	100 pcs.
502-156	SAAR2156	Fluorhibit	454g
502-167	SA84994401	Tin Capsules, 8mm, 8x20mm	100 pcs.
502-173	SAAR173	Alphacell II - HP	750g
502-174	SAAR2174	Sansorb II, CO2 Absorb, 20-30 Mesh	500g
502-176	SAAR2176	Sansorb II, CO2 Absorb, 8-20 Mesh	500g
502-177	SAAR2177	Fine Quartz Wool	50g
502-181	SA999959	Tungstic Oxide, Granular	60g
502-183	SAAR2183	Nickel Capsule, 6.6 x 10.0mm	100 pcs.
502-184	SA990747	Benzoic Acid Powder, 99%	1g
502-186	SA76152301	Tin Foil Cups	100 pcs.
502-187	SA84992101	Tin Capsule, 3.5mm x 9mm	100 pcs.
502-188	SA502188	Alumina Oxide Pellets, FP428	200g
502-189	SA99060302B	Copper Sticks	100g
502-190	SA99060301	Copper Oxide, Fine Wires	100g
502-203	SAAEB2031	Stearic Acid Standard	1g
502-204	SAAEB2034	Dinitrobenzoic Acid Standard	1g
502-205	SA999946	Caffeine Standard	1g
502-206	SA502206W	Silver Capsules, 3.2 x 4 mm	1000 pcs.
502-207	SA33840018	Cystine Standard	1g
502-208	SA990750	Sulphanilamide Standard	g1
502-209	SAAR2209	Sulphamethazine Standard	5g
502-211	SAAR211	Glycine 99%	50g
502-213	SAAR2213	Leak Detection Solution	113g
502-225	SA990797	Silver Tungstate/MGO	50g
502-226	SA990789	Silver Vanadate	15g
502-227	SA502227W	Tin Capsules, 3.2mm x 4mm	500 pcs.
502-231	SAAR673	Wrought Iron Chip Standard 0.0004% C, 0.0017% S	454g
502-258	SA501621W	Copper Turnings, Degassed	200g
502-266	SAAR205	Firing Fuse Wire, 300 Determinations X 0.006" Dia	
502-272	SAAR2017	Corn Gluten, C,H,N,S	30g
502-273	SAAR2018	Alfalfa, C,H,N,S	30g
502-274	SAAR2019	Wheat Flour, C,H,N,S	30g
502-275	SAAR2020	Rye Flour, C,H,N,S	30g
502-276	SAAR2026	Oat Meal, C,H,N,S	30g
502-277	SAAR2027	Barley, C,H,N,S	30g
502-278	SAAR2028	Rice Flour, C,H,N,S	30g
502-282	SA99060302B	Copper Sticks	100g
502-295	SA501621W	Copper Turnings, Degassed	200g
502-297	SAAR008B	Alphacel II, HP	13.6kg
502-298	SAAR298	Sulfamethazine	25g
502-301	SA999959	Tungstic Oxide, Granular	60g
502-304	SA502304E	Copper Sticks in Argon Vial	100g
502-304-500	SA502304500E	Copper Sticks in Argon Vial	500g
502-308	SAAR4016	Soil Standard, 2% C & S	25g
502-310	SAAR2310	Steel Wool Filter Material	454g
502-311	SAAR2311	Furnace Reagent 1000	1 pc.
502-318	SAAR4017	Soil Standard, C=0.5%, S=0.5%	
502-319	SAAR4005	Ore Standard, 30g, C=1.42%, S=1.41%	
502-320	SAAR4006	Ore Standard, 30g, C=4.40%, S=4.24%	
502-321	SAAR321	Alphacat Accelerator	250g
502-335	SAAR2335	Foil Crucible Liner, TGA	250 pcs.
502-338	SAAR338	Alphagel Capsules, for FP428	400 pcs.
502-343	SA2343	Nickel Boat Liner,	100 pcs.
502-344	SA502344E	Nickel Baskets, Low Oxygen, 1g	100 pcs.

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
502-345	SAAR2345	Nickel Baskets, Low Oxygen, 1.5g	100 pcs.
502-351	SAAR351	Halogen Scrubber	50g
502-359	SAAR2359	Furnace Reagent, MGO	200g
502-372	SAAR4007	Ore Standard, 30g, C=7.62%, S=3.25%	
502-374	SAAR2374	Heater Particulate Filter	50 pcs.
502-382	SAAR339	Medium Alphagel Caps	400 pcs.
502-397	SAAR2397	Large Tin Foil Cups	100 pcs.
502-397-400	SA502397400	Large Tin Foil Cups	400 pcs.
502-435	SAAR1704	Coal Standard, 1% Sulfur	
502-489	SAAR2489	Porous Crucible, CHN	10 pcs.
502-493	SAAR2156	Fluorhibit	454g
503-007	SAAR007	Detachable Float Valve	
503-020	SAAR3020	Boron Distillation Apparatus	
503-031	SAAR905	Funnel (3 Leg)	
503-032	SAAR901	Scoop	
503-033	SAAR033	Sulfur Trap, Filled	
503-051	SAAR051	Oxygen Baffle for 1 1/8" ID Tube	
503-052	SAAR052	Oxygen Baffle for 1 1/4" ID Tube	
507-010	SAAR911H	Catalyst Tube (Hopcolite)	
507-011	SAAR911E	Catalyst Tube (Unfilled)	
507-012	SAAR012	Reducer	
507-100	SAAR7100	Catalyst Furnace, Rebuilt	
507-200	SAAR7200	Catalyst Furnace, Rebuilt	
507-207	SAAR911	Catalyst Tube, Filled, Platinized Silica	
507-207	SAAR911R	Catalyst Tube, Filled, Rare Earth	
508-300	SAAR8300	Catalyst Heater, Rebuilt	
508-400	SAAR8400	Catalyst Heater, Rebuilt	
508-500	SAAR8500	Catalyst Heater, Rebuilt	
508-600	SAAR8600	Catalyst Heater, Rebuilt	
511-036	SAAR182	4mm ID Evacuated Pin Tube	
511-038	SAAR1038	Pin Tube Holder	
511-039	SAAR181	3mm ID Evacuated Pin Tube	
511-390	SAAR183	6mm ID Evacuated Pin Tube	10 pcs.
516-000	SAAR906	Purifying Train	
516-004	SAAR504	Rotometer	
516-005	SAAR505	Drying Tower	
516-006	SAAR506	Acid Tower	
516-063	SAAR507	Needle Valve Complete	
517-204	SAAR204	Rubber Squeeze Bulb	
518-016	SAAR016	Burette	
518-024	SAAR024	Light Director	
518-030	SAAR017	Vessel	
519-005	SAAR900	Ignitor	
519-009	SAAR009	Pedestal	
519-044	SAAR004	Hydrocarbon Sulfur Tube	
521-112	SAAR114	Silicone Tubing, 4" x 3/8" ID	
521-113	SAAR113	Silicone Ring Seal	
521-186	SAAR187	Solid State Rectifier	
521-287	SAAR287	Suppressor Coil	
521-305	SAAR932	Drying Tube	
528-018	SAAR3818	Ceramic Crucible	
528-028	SAAR8028	Filtering Crucibles, Small	
528-030	SAAR8030	Filtering Crucibles, Large	
528-040	SAAR040	Cover, 4mm Hole	
528-042	SAAR042	Porous Covers	
528-043	SAAR9880	Cover, 10mm Hole	1000 pcs.

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
528-050	SAAR3818B	Ceramic Crucibles, 50/bag	500 pcs.
528-053	SAAR8053	Combustion Boat, 4" x 1/2"	
528-120	SAAR8120	Low Sulfur Crucible	
528-203	SAAR8203	Combustion Boat For SC132-32-CR-12	500 pcs.
528-203-250	SAAR8203250	Combustion Boat for SC132/32/CR-12	250 pcs.
528-230	SAAR5282	Combustion Boat, Glazed for Liquids	50 pcs.
528-313	SAAR8313	Combustion Tube 1 1/8" x 30" OBE	
528-314	SAAR8314	Combustion Tube 1 1/4" x 30" OBE	
528-331	SAAR8331	Combustion Tube 1 1/8" x 27" TE	
528-332	SAAR8332	Combustion Tube 1 1/8 x 30" TE	
528-333	SAAR8333	Combustion Tube 1 1/4" x 30" TE	
529-007	SAAR529	Ceramic Pedestal, 3.05"	
529-011	SAAR9011	Ceramic Liner SC32/132/CR-12	
529-012	SAAR9012	Boat Stop Kit (2 x AR9012)	
529-018	SAAR3203	Ceramic Nozzle 320-00203 Lecomelt	
529-022	SAAR9022	Boat Stop Half, Honeycomb Fine	
529-026	SAAR026	Crucible CHN-600	50 pcs.
529-028	SAAR9023	Boat Stop Half SC132/32	
529-031	SAAR9031	Cone Support D	
529-042	SA529042W	Ceramic Crucible, MAC-400/TGA, Small	5 pcs.
529-043	SA529043W	Ceramic Cover, MAC-400/TGA, Small	5 pcs.
529-046	SAAR9046	Ceramic Tube TF-10, 1 1/4" x 23" L	
529-047	SA529047W	Ceramic Crucible, MAC-400 Large	5 pcs.
529-048	SA529048W	Ceramic Cover, MAC-400 Large	5 pcs.
529-054	SAAR9023	Boat Stop Half SC32/132	
529-055	SAAR026	Crucible CHN	50 pcs.
529-203	SAAR8203P	Combustion Boat SC32/132/CR12	
529-204	SAAR9204P	Combustion Boat, S-144DR	50 pcs.
529-204-250	SAAR9204250	Combustion Boat, 250/box	250 pcs.
532-065	SAAR065	Burette for 532-500	
532-066	SAAR066	Vessel for 532-500	
532-500	SAAR532	Automatic Sulfur Titrator	
533-104	SAAR801	Starch Dispensing Pipette	
537-103	SAAR103	1/4" OD Polyflow Tubing	
540-140	SAAR140	Bulb	
544-147	SAAR4147	Leveling Bottle	
549-007	SAAR907	Armature	
549-072	SAAR9072	Restrictor	
549-076	SAAR176	Bulb, Light Source for 518, 532	
550-004	SAAR404	Oscillator Tube, EIMAC	
550-044	SAAR044	Insulator, Porcelain 521/523	
550-046	SAAR046	Resistor, 521/523	
550-047	SAAR0047	RF Choke for 521/523	
550-053	SAAR0053	Power Switch for 521/523	
550-119	SAAR019	Right Angle Bend	
550-120	SA550120E	Sulfur Combustion Tube	1 pc.
550-122	SA550122E	Carbon Combustion Tube	1 pc.
550-123	SAAR123	Pedestal Insert	
550-124	SAAR124	Pedestal, 2.50"	
550-165	SAAR165	Spring Clip for Combustion Tube	
550-167	SAAR550	Furnace Fan for Models 521/523, 115VAC	
550-182	SAAR180	QEGC Crucible	
550-183	SAAR179	Quartz Crucible Insert for AR180	
550-187	SAAR186	Capacitor	
566-008	SAAR006	Induction Coil	
568-016	SAAR428	"O" Ring	1 pc.

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
572-016	SAAR72016	Control Valve for 572	
572-029	SAAR029	Float Valve Assembly	
572-031	SAAR2031	Splashproof Device	
578-672	SAAR672	"O" Ring	1 pc.
585-235	SAAR200	AL203 Wheel, 7" x 3/4" H x 1/16"	
585-270	SAAR270	AL203 Cut-off Wheel 7 x 1" x 1/16"	
585-359	SAAR5359	Pin Cutter Spacer, for 1/2g Pin	
585-360	SAAR5360	Pin Cutter Spacer, for 1g Pin	
585-500	SAAR2000X	Pin Cutter, 1 H.P., 230VAC, 50 Hz	
585-600	SAAR2000	Pin Cutter, 1 H.P., 115VAC, 60 Hz	
589-305	SAAR9305	1/4" Rubber Sleeve	
589-364	SAAR9364	Oven Fan, WR12, TC30, Dayton	
589-551	SAAR422	"O" Ring (Auto Drop)	1 pc.
589-859	SAAR859	Thermistors, WR-12, 750	
589-887	SAAR887	"O" Ring	1 pc.
589-891	SAAR9891	Timer Relay 750/589	
593-100	SAAR593	Timer, Sulfur Titrator	
598-142	SAAR859	Thermistors, WR-12, 750	
598-157	SAAR157	Moisture Trap, Filled, 750, 598	
598-158	SAAR158	Moisture Trap, Unfilled, 750, 598	
598-257	SAAR257	Replaced by 614-591	
598-614	SAAR613	Bulbs - Switches 6.3V (750-100)	
601-064	SAAR064	Quartz Heater Tube, CF-10	
601-065	SAAR1065	Pyrex Reagent Tube, CF-10	
601-132	SAAR132	Combustion Tube CHN 800	
601-133	SAAR1133	"O" Ring (1.0 x 1.25) CHN	1 pc.
601-158	SAAR1158	"O" Ring CHN 800	1 pc.
601-159	SAAR1159	"O" Ring CHN 800	1 pc.
601-201	SAAR1201	Reaction Tube DH103	
601-390	SAAR390	Quick Disc Tube FP428	
601-397	SAAR397	Shim, Combustion Tube FP428	6 pcs.
601-437	SA601437E	Quartz Catalyst Tube, FP428	
601-438	SA990759	Combustion Tube, for NP228/428, SN 1 thru 3445	
601-441	SAAR441	Quartz Lance Tube, FP428	
601-442	SAAR442	"O" Ring, FP428, RC612	1 pc.
601-482	SAAR482	Honeycomb Support, FP428	
601-483	SAAR483	Printer Paper, 9 1/2 x 11	600 pcs.
601-484	SAAR484	Printer Ribbon Cartridge	
601-490	SAAR1490	Pelletizing Foiler Device	
601-504	SAAR1504	"O" Ring, Loading Block, FP228	1 pc.
601-505	SAAR1505	"O" Ring, Combustion Tube, FP228, RC612	1 pc.
601-506	SA601506E	Combustion Tube, FP228	
601-684	SAAR684	"O" Ring, Loading Head Assembly TC600	1 pc.
601-691	SAAR1691	.125 Cflex Tubing, 12"	
601-691-030	SAAR169130	Cflex Tube, 30"	
601-691-036	SAAR169136	Cflex Tube, 36"	
601-862-110	SAAR1862	Furnace, FP-428, 550-650 Watts	
601-920	SAAR920	"O" Ring	1 pc.
601-939	SAAR1939	"O" Ring, Doser, FP428	1 pc.
601-981	SAAR981	Side Arm Cleaning Brush, FP428	
602-615	SAAR2615	Furnace Housing w/o End Caps FP2000	
602-802	SAAR802	"O" Ring, RH 402, RC612	1 pc.
602-803	SAAR803	"O" Ring, RH402	1 pc.
602-969	SAAR969	5000 Determinatons Package for FP428	
603-047	SAAR63047	Ureth Tubing, .061 x .12 OD, Red, Per Foot	
603-048	SAAR63048	Ureth Tubing, .061 x .12 OD, Blue, Per Foot	

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
603-049	SAAR63049	Ureth Tubing, .061 x .12 OD, Green, Per Foot	
603-050	SAAR63050	Ureth Tubing, .061 x .12 OD, Black, Per Foot	
603-051	SAAR63051	Ureth Tubing, .061 x .12 OD, White, Per Foot	
603-052	SAAR63052	Ureth Tubing, .061 x .12 OD, Orange, Per Foot	
603-493	SAAR3493	"O" Ring	1 pc.
603-706	SAAR706	Pedestal for RH402	
603-707	SAAR707	Reaction Tube RH402	
603-710	SA603710E	Quartz Crucibles, RH402	24 pcs.
603-910	SAAR3910	"O" Ring	1 pc.
603-970	SAAR3970	"O" Ring RF402, Upper Load Mech	1 pc.
603-971	SAAR3971	"O" Ring RH402, Upper Load Mech	1 pc.
604-112	SAAR604	Quartz Catalyst Tube	
604-114	SAAR923	Pedestal Only, Quartz with Fitting	
604-163	SAAR163	Brass Brush	
604-172	SAAR172	"O" Ring	1 pc.
604-241	SAAR4241	Reaction Tube, HF402/RH402	
604-271	SAAR271	SS Fluted Brush, Dia Slotted	
604-373	SAAR4373	Capsule Holder	
604-496	SAAR4496	Teflon Flx Tubing,.125 x .032W, ft	
604-518	SAAR518	SS Cleaning Brush, Unslotted	
604-541	SAAR4541	Lance Tube CHN 900/932	
605-520-087	SAAR6055	"O" Ring Package, CS200/300/400	1 pc.
605-814	SAAR5814	Pedestal Adapter Assembly	
605-815	SAAR5815	Pedestal Adapter, CS144-CS444	
605-816	SAAR816	Ceramic Pedestal, CS144-444 SiN	
605-817	SAAR817	Pedestal Adapter, CS125	
605-818	SAAR818	Pedestal Adapter 606/816, EC12, HF10/20	
605-828	SAAR828	"O" Ring, FP428, SN 1 thru 3445	1 pc.
605-830	SAAR830	"O" Ring CHN	1 pc.
605-875	SA605875	Furnace Red/Ox Tube, CHN900	
606-202	SAAR202	Furnace Fan for FP2000 Series	
606-224	SAAR6224	Thermocouple, SC432 Series	
606-286	SAAR563	Short Strap, SC432	
606-308	SAAR6308	Stepped Ceramic Boat Stop, SC432	
606-309	SA606309E	Outer Combustion Tube SC432	
606-310	SAAR1477	"O" Ring Combustion Tube, SC432	1 pc.
606-312	SAAR1476	"O" Ring Combustion Tube, SC432	1 pc.
606-313	SA606313E	Inner Combustion Tube, SC432	
606-316	SAAR6316	Lance Tube Assembly, SC-432	
606-317	SAAR6317	Stainless Purge Tube	
606-318	SAAR6318	Heating Elements (3), SC432	
606-319	SAAR6319	Lead Wire Braided 6.5, SC-432	
606-320	SAAR6320	Braid Clamp for SC432 Element	
606-327	SAAR6327	Gas Drying Tube SC432	
606-333	SAAR6333	"O" Ring, Lance, SC432	1 pc.
606-340	SAAR6340	Drilled Lance Tube, FP2000	
606-403	SAAR403	Collection Tube Brush, SC432	
606-444	SAAR6644	#2 Stepped Ceramic Boat Stop	
606-547	SAAR6547	Boat Puller, SC432 Autoloader	
606-550	SAAR6550	Finger Arm for AR6547	
606-601	SAAR6601	Heat Elements (4) CS432, CS144	
606-601-100	SAAR6601100	Heating Element Kit	
606-602	SAAR6602	Furnace End Cap FP2000	
606-655	SAAR6655	Purifier Tube, SC432	
606-747	SAAR6747	2 Hole Lance Tube, FP2000, 245mm	
608-339	SAAR8339	Combustion tube, Old Syle, CHN1000, 14mm OD x 200mm	

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
608-342	SAAR342	Combustion Tube, Quartz, CH1000	
608-343	SAAR8343	Furnace Element, CHN1000	
608-344	SAAR344	Combustion Tube, Outer, CHN1000	
608-378	SAAR1476	"O" Ring, CHN1000	1 pc.
608-379	SAAR8379	Quartz Wool Strips, 15"	10 pcs.
608-505	SAAR8505	Combustion Tube, Ceramic, CHN1000	
608-866	SAAR866	Cflex Tubing .187"/foot	
609-114	SAAR9114	"O" Ring, Side Arm, SN3446 up	1 pc.
609-116	SA609116E	Combustion Tube, FP428, SN3446 up	
609-129	SA609129E	Lance, Straight	
609-170	SAAR9170	Quick Cap Kit	
609-171	SAAR9171	Plunger for Cellulosic Tubes	
609-172	SAAR9172	Tubing Cellulosic, .210ID	20 pcs.
609-271	SAAR9271	Radial Shaft Seal, FP2000 Autoloader	
609-394	SAAR394	Heat Shield for FP/NS2000	
609-474	SAAR9474	Finger Arm for AR1207 Shaft	
609-728	SAAR9728	Shaft Sleeve Bearing, FP2000 Autoloader	
610-735	SAAR735	Furnace Brush, CS444, Dia Slotted	
611-207	SAAR1207	Puller Shaft Assembly with Finger	
611-255	SAAR9171	Plunger Assembly	
611-351-189	SAAR691	Electrode Insert	
611-351-223	SAAR1351	Combustion Tube, HD	
611-431	SAAR1431	"O" Ring 132	1 pc.
611-476	SAAR1476	"O" Ring, Inner Combustion Tube	1 pc.
611-477	SAAR1477	"O" Ring, Outer Combustion Tube	1 pc.
611-790	SAAR790	"O" Ring	1 pc.
612-340	SAAR2340	"O" Ring, Thermistor Block, Truspec	1 pc.
613-563	SAAR563	Short Strap, SC432	
614-591	SAAR257	Cabinet Fan, 110V	
614-822-100	SA874000181W	Nickel Boat, Small for AMA 254 Hg	1 pc.
614-822-101	SA874000180W	Small Boat Loop	1 pc.
614-822-114	SA874001181W	Nickel Boat, Large for AMA 254 Hg	1 pc.
614-822-133	SA874001180W	Large Boat Loop	1 pc.
614-961	SAAR4961	Large Crucible,FP528 Truspec	10 pcs.
614-962	SAAR4962	Lance, FP528	
616-063	SAAR6063	Boat Puller,Sulfur/Carbon	
616-068	SAAR6068	"O" Ring, Inner Combustion Tue S/SC144	1 pc.
616-070	SAAR6070	Inner Combustion Tube S/SC-144	1 pc.
616-086	SAAR6086	Stepped Boat Stop, S/SC-144	
616-089	SAAR6089	Outer Combustion Tube, S/SC-144	1 pc.
616-108	SA616108	Combustion Tube, FP528	1 pc.
616-138	SAAR6138	"O" Ring, FP-528	1 pc.
616-139	SAAR6139	"O" Ring, FP-528	1 pc.
616-146	SAAR6146	Lance, FP528	
616-504	SAAR6504	Particle Filter, FP528	
616-605	SAAR7605	Large Crucible for Liquid Samples	10 pcs.
616-622	SAAR6622	Micro Switch, CN 2000; Switch Only No Leads	
616-632	SAAR6632	Lance, SC-432 & FP-2000, 253mm	
616-635	SAAR394	Heat Shield for FP/NS2000, Replaces 609-394	
617-083	SAAR7083	"O" Ring, in Scrubber Cat Heater	1 pc.
617-136	SAAR7136	"O" Ring, Loading Head, Truspec	1 pc.
617-440	SAAR440	"O" Ring, Filter Assembly TC600	1 pc.
617-441	SAAR7441	Micro Screen Filter, TCH600	
617-548	SAAR7548	Filter Tube, TC500, RC612	
617-716	SAAR7716	Quartz Filter Tube, TC500	
618-376	SAAR376	Electrode Insert	

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
618-482	SAAR691	Electrode Insert	
618-603	SAAR603	Catalyst Tube, 13.0"	
618-657	SAAR6186	FTG Adapter, RC612	
618-809	SAAR8809	Reagent Tube, TC600	
619-065	SA619065E	U-Shaped Combustion Tube for Truspec	1 pc.
619-154	SAAR9154	Catalyst Heater Tube for Truspec	1 pc.
619-192	SAAR9192	"O" Ring, Bottom Load Head for Truspec	1 pc.
619-233	SAAR9233	Thermocouple for Truspec	
619-268	SAAR9268	High Sensitivity Reagent Tube for Truspec	
619-316	SAAR9316	Large Volume Reagent Tube for Truspec	
619-425	SAAR9425	Straight Reagent Tube or Truspec	
619-455	SAAR9455	Large Reagent Tube	
619-566	SAAR9566	Braided Strap	
619-591-699	SAAR1699	Particle Filter	1 pc.
619-774	SAAR9774	Lance, Truspec, CHNS Micro	
619-776	SAAR9776	Reagent Tube 233mm	
619-777	SAAR9777	Directive Middle Tube CHNS	
619-803	SAAR3302	Upper Electrode Brush, 21 mm diameter	
619-880	SAAR9880	10mm Hole Cover	1000 pcs.
619-895	SAAR619	Graphite Crucible, RHN600	1000 pcs.
619-895-100	SAAR619100	Graphite Crucible, RHN600	100 pcs.
620-814	SAAR691	Electrode Insert	
621-331	SAAR9047	TGA Crucible	5 pcs.
625-067	SAAR5067	Elements	6 pcs.
625-068	SAAR25068	Element Clips	
625-501	SAAR6255	Thermocouple, RC612	
625-505-125	SAAR5125	Catalyst Heater Tube RC612	
625-505-138	SAAR5138	Heating Element Assembly RC612	
625-505-168	SAAR5168	Combustion Tube, RC612	
625-505-271	SAAR5271	Knurled Adapter, RC612	
625-510-402	SAAR5104	Combustion Tube, RC612	
625-511-154	SAAR511154	Lower Crucible, CHN	5 pcs.
660-010-130	SAAR9233	Thermocouple for Truspec Type K	
701-030	SAAPT1030	Torch Assembly, ICP-2000 LECO	
702-557	SAAR2557	"O" Ring, Truspec	1 pc.
750-301	SAAR933	Drying Tube, 9.0"	
751-300	SAAR1300	Electronic Balance, 4th Place	
760-034	SAAR034	Graphite Crucibles	
760-072	SAAR072	28 DVC Bulbs	
760-138	SAAR138	Tweezers, Straight	
760-139	SAAR766	Tweezers, Curved Deluxe	
760-144	SAAR022B	Brush, 90" Bend	
760-145	SAAR0145	Powder Funnel for Graphite Crucibles	
760-193	SAAR193	Metal Screen	
760-224	SAAR224	"O" Ring	1 pc.
760-225	SAAR425	"O" Ring (Loading Head Plunger)	1 pc.
760-226	SAAR428	"O" Ring	1 pc.
760-293	SAAR293	Brush, Cleaning, Small	
760-304	SAAR0304	Teflon Washer	
760-316	SAAR316	Washer	
760-317	SAAR317	Glass Tube, 6.64"	
760-333	SAAR333	Drying Tube, 8"	
760-335	SAAR424	"O" Ring	1 pc.
760-359	SAAR359	"O" Ring (Quick Disconnect)	1 pc.
760-414	SAAR414	Graphite Crucibles	1000 pcs.
760-841	SAAR194	"O" Ring	1 pc.

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
761-047	SAAR047	3" Drying Tube	
761-133	SAAR133	Drying Tube, 17.125"	
761-190	SAAR190	Drying Tube, 11.75"	
761-383	SAAR293H	Small Brush Handle Only	
761-492	SAAR492	"O" Ring FP428	1 pc.
761-618	SAAR618	Trap Repair Charge/WR-12	
761-651	SAAR1609	"O" Ring, Lance, FP/CHN	1 pc.
761-739	SAAR639	Tin Flux	500g
761-747	SAAR611	Schutze's Reagent	100g
761-918	SAAR918	"O" Ring	1 pc.
761-929	SAAR929	Crucible Tongs for AR3818	
762-023	SAAR023	Nut Adaptor Seal for CS44/IR12	
762-056	SAAR056	Combustion Tube	1 pc.
762-057	SAAR416	"O" Ring (Dust Filter)	1 pc.
762-058	SAAR058	"O" Ring	1 pc.
762-059	SAAR614	"O" Ring (Cup Assembly)	1 pc.
762-060	SAAR417	"O" Ring Backup	1 pc.
762-062	SAAR913	Oxygen Lance	
762-192	SAAR192	IR So2 Cell, CS44/IR32/33	
762-357	SAAR357	"O" Ring (Upper Adapter)	1 pc.
762-501	SAAR2501	Handle for Rebuilt 521	
762-523	SAAR523	Oven Fan Replaces Howard Fan	
762-621	SAAR2621	IR Co2 Cell, CS44/IR12	
762-627	SAAR2627	Biscuit Fan, CS44	
762-695	SAAR695	Tin Metal	
762-697	SAAR697	Large Throat Brush	
762-741	SAAR650	0.1g H/O/N in Titanium Standard 0.0061% H, 0.1555% O, 0.0101% N	100 pcs.
762-747	SAAR546	1g Hydrogen in Steel Pin Standard 1.1 PPM H	100 pcs.
762-823	SAAR697H	Brush with Handle	
762-936	SAAR936	Glass Scoop for V205, Extra Large	
763-026	SAAR027B	Alphacel III Tungsten 30# - 20 + 40	
763-027	SAAR266B	Alphacel Tungsten 30# - 12 + 20	
763-028-500	SAAR027V	Alphacel III, Tungsten	227.27kg
763-029	SA711526	Nickel Baskets 1.5g	100 pcs.
763-065	SA711525	Nickel Baskets 1g	100 pcs.
763-191	SAAR103	Tubing Polyflo, 1/4" OD	
763-213	SAAR213	Graphite Crucible for RH-1E	
763-261	SAAR261	Filter Paper	
763-263	SAAR266	Alphacel 5# - 12 + 20	
763-265	SAAR265	Glass Wool	227g
763-266	SAAR027	Alphacel Tungsten, - 20 + 40	2.3kg
763-295	SAAR295	Thermostat, Upper Electrode	
763-379	SAAR3379	"O" Ring, FP428	1 pc.
763-467	SAAR467	Iron Chip 25#	11.4kg
763-873	SAAR3873	Isolation Box Bulb, 28VDC	
764-015	SAAR015	Bulb for 532-500, 28VDC	
764-078	SAAR4078	1g Calibration Weight	
764-098	SAAR098	Cellulose Material	454g
764-216	SAAR4216	Helium Regulator Assembly	
764-242	SA764242	Tin Pellets	454g
764-490	SAAR490	Pedestal, 2-3/4"	
764-502	SAAR4544	Rotron Furnace Fan	
764-545	SAAR2832	Lube Oil, 4 oz, 0.25% Sulfur	
764-546	SAAR2835	Lube Oil, 4 oz 1.00% Sulfur	
764-547	SAAR2836	Lube Oil, 4 oz 2.00% Sulfur	

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
764-640	SAAR334	100 PF Capacitor	
764-641	SAAR335	200 PF Capacitor	
764-720	SAAR064	Quartz Disconnect Tube, CF-10	
764-762	SAAR762	Glass Tube 11.88" NP28	
764-775	SAAR775	WR-12 Trap Material	
765-068	SAAR5068	IR Oxygen Cell for RO17	
765-523	SAAR523	Oven Fan, Howard Kit	
765-717	SAAR917	Oxygen Lance (WR 12)	
765-746	SAAR746	Rotometer, WR12	
765-818	SAAR5818	Optic Isolation Cable	
765-933	SAAR5933	Lubricool Cooling Agent	
765-976	SAAR5976	"O" Ring, AC AC100, FP428	1 pc.
766-036	SAAR6036	Pneumatic (Air) Regulator Assembly	
766-053	SAAR766	Deluxe Crucible Tweezers, Curved	
766-557	SAAR6557	"O" Ring	1 pc.
766-646	SAAR6646	Air Filter Coating	
766-823	SAAR823	Male Cliplock Fitting	
767-243	SAAR243	Female Cliplock Fitting	
767-277	SAAR277	Graphite Crucibles	
767-385	SAAR085	Quartz Wool, Coarse	28g
767-500	SAAR914	Lance, IR33	
767-518	SAAR7518	Retaining Wire Clip Spring	
767-541	SAAR541	Drying Tube	
767-542	SAAR542	Scoop for Tungsten, 0.01cc, Small	
767-757	SAAR1706	Coal Standard, 2.0% Sulfur	
767-758	SAAR1707	Coal Standard, 2.5% Sulfur	
767-759	SAAR720	Coke Standard, 0.85% Sulfur	
767-838	SAAR838	"O" Ring	1 pc.
768-774	SAAR9880	10mm Hole Cover	1000 pcs.
768-798	SAAR8798	Clean Out Plug for SC32/132/CR	
768-944	SAAR944	"O" Ring	1 pc.
768-980	SAAR980	Particle Filter - SC32/132/RC612	
769-335	SAAR469	SS Relay, HF10, 44F 6848 240 D45	
769-520	SAAR9520	Graphite Crucibles/RH-2	
769-608	SAAR608	Halogen Trap Material	200g
769-610	SAAR610	Halogen Trap Material	225g
769-761	SAAR761	Crucibles, RH-3	100 pcs.
769-821	SA990713	Molecular Sieve for RH-2	
769-898	SAAR898	Lower Pedestal	
769-994	SAAR414	Graphite Crucible	1000 pcs.
770-247	SAAR247	Filter for SC32/132	4 pcs.
770-436	SAAR523	Oven Fan with Adapter	
770-764	SAAR764	Printer Paper	
770-778	SAAR778	Ribbon Spool Assembly	
771-015	SAAR1015	Manual WT Comp Control CS46	
771-030	SAAR030	Bubble Tube	
771-040	SAAR3913	"O" Ring	1 pc.
771-275	SAAR1275	Pneumatic Cylinder, CS46/IR112	
771-354	SAAR354	Oscillator Tube CS-144-444	
771-367	SAAR1367	Quartz Crucibles, DC12	50 pcs.
771-373	SAAR373	Drying Tube	
771-401	SAAR401	Reaction Tube for RH-2	
771-407	SAAR407	"O" Ring	1 pc.
771-425	SAAR1425	Power Triode RH-2, 7.5V	
771-446	SAAR429	"O" Ring	1 pc.
771-606	SAAR606	Pedestal CS-46 & EC-12	

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
771-883	SAAR083	Bubble Tube (EC 12)	
771-991	SAAR991	"O" Ring 028	1 pc.
772-005	SAAR760	Chart Paper Roll, RH2/3	
772-269	SAAR269	"O" Ring	1 pc.
772-381	SAAR2381	Lance Cleaning Wire	
772-402	SAAR2402	Lance SC-32/132/CR-12	
772-403	SAAR6403	Drilled Lance Tube, SC432	
772-446	SAAR446	Glass Viewing Shield, SC32/132	
772-449	SAAR449	Rod, Extensions, Glass Shield	
772-460	SAAR2460	Pressure Sensor	
772-469	SAAR469	Replaced by 781-165 (AR469)	
772-508	SAAR962	Element for SC32/132/CR12	
772-520	SAAR520	"O" Ring SC32, CS244, 344, RC612	1 pc.
772-522	SAAR7073	Thermocouple Only	
772-560	SAAR2560	Chopper Motor	
772-568	SAAR233Q	Ceramic Adhesive Alumina	
772-738	SAAR738	"O" Ring	1 pc.
772-739	SAAR544	Catalyst Tue (10.75")	
772-747	SAAR2747	Fan for CS244	
772-795	SAAR795	Oven Fan, CS244 & TC136	
772-905	SAAR057	Combustion Tube, CS 144-444	
772-906	SAAR2906	Combustion Tube Ring Holder	
772-910	SAAR2910	"O" Ring, Piston Seal, FP228	1 pc.
773-106	SAAR106	Delivery Tube	
773-136	SAAR3136	Fitting, SS for FP428	
773-141	SAAR320	"O" Ring	1 pc.
773-311	SAAR311	SST Filter	
773-314	SAAR3314	Unfluted Cleaning Brush	
773-316	SAAR789	Lance Tube, CS125	
773-320	SAAR320	"O" Ring	1 pc.
773-561	SAAR561	SS Cleaning Brush, Straight Slotted	
773-579	SAAR579	Metal Scoop	
773-662	SAAR3662	Polyflo Flex Tubing, CHN	
773-913	SAAR3913	"O" Ring	1 pc.
773-928	SAAR3928	Delivery Tube for SC32/132	
773-931	SAAR9931	Teflon Disconnect Block, SC132	
773-932	SAAR3932	Combustion Tube Repair Service	
774-031	SAAR4031	Quick Disconnect Tube 1.38" OD	
774-174	SAAR4174	Purge Tube	
774-204	SAAR4204	AC Crucibles, SST	10 pcs.
774-205	SAAR205	Firing Wire, 300 Determinations x 0.006" Dia	
774-208	SAAR208	Benzoic Acid Tablet	10 pcs.
774-211	SAAR4211	"O" Ring 2.6 x 2.875	1 pc.
774-253	SAAR4253	"O" Ring	1 pc.
774-283	SAAR4283	"O" Ring	1 pc.
774-425	SAAR432	"O" Ring, Closure Block	1 pc.
774-460	SAAR460	Thermocouple Assembly SC32/132/CR12	
774-489	SAAR489	"O" Ring, Upper Electrode	1 pc.
774-570	SAAR4570	Vacuum Pump for Early SC32 110VAC	
774-571	SAAR4571	Vacuum Pump for SC32/132, 24VDC	
774-613	SAAR4613	"O" Ring, HF400 Pinch Valve	1 pc.
774-638	SAAR638	Disco Tube, 6" Straight, R0116	
774-649	SAAR653	Cleaning Brush with Handle	
774-652	SAAR652	Brush for TC-136	
774-690	SAAR690	Knurled ring Collar, EF100/EF500	
774-691	SAAR691	Replaced by 618-482	

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
774-692	SAAR692	Lower Electrode Exchange/Repair Only	
774-809	SAAR809	Induction Coil HF10, S/Oil	
774-841	SAAR841	Thermistor Sensor, Lower Drain	
775-288	SAAR3403	Benzoic Acid Powder	454g
775-306	SAAR5306	Secondary Filter	
775-312	SAAR312	"O" Ring, Lower Electrode	1 pc.
775-313	SAAR313	Teflon Backup Ring	
775-314	SAAR314	"O" Ring, Piston Block	1 pc.
775-330	SAAR5330	Four Way Valve	
775-355	SAAR355	Catalyst Tube, 10.75"	
775-416	SAAR293	Cleaning Brush, TC-136, TN-14	
775-416	SAAR293H	Brush Handle	
775-431	SAAR5431	Inner Graphite Crucible	1000 pcs.
775-433	SAAR433	Outer Graphite Crucible	50 pcs.
775-446	SAAR3930	Electrode Insert Repair Service	
775-601	SAAR601	Bubble Tube, Left	
775-625	SAAR625	Momentary Push Switch	
775-677	SAAR677	"O" Ring	1 pc.
775-814	SAAR402	Reaction Tube, Thick Wall, RH2/3	
775-859	SAAR5859	Wire Cleaning Brush, SC32/132	
775-880	SAAR880	Quick Disconnect Tube, 4"x, Straight,FP228	
775-892	SAAR431P	Inner Graphite Crucible	100 pcs.
776-017	SAAR6017	"O" Ring, Ballast Plate, Truspec	1 pc.
776-095	SAAR095	Silicone Tubing Dust Coll.	
776-129	SAAR129	Induction Coil, CS144-444	
776-234	SAAR6234	10 Micron Filter, 1gpm	
776-247	SAAR053	Graphite Crucibles - Discontinued!!! Replaced by SAAR6247	
776-247	SAAR6247	Graphite Crucibles	1000 pcs.
776-265	SAAR6265	Adjustable Pedestal Assembly, 923/4060/6922	
776-279	SAAR279	Catalyst Heater Tube	1 pc.
776-283	SAAR283	Quartz Lance CHN-600/800	1 pc.
776-284	SAAR284	"O" Ring, Combustion Tube FP228	1 pc.
776-286	SAAR286	Combustion Tube CHN 600	1 pc.
776-685	SAAR6685	Tyna-Mite Solenoid	
776-923	SAAR923	Quartz Pedestal Post CS144-444	
776-931	SAAR606	Pedestal HF10/20	
776-962	SAAR962	Element for SC32/132/CR12	
776-978	SAAR978	Glass Scoop, MAC 400	
776-992	SAAR6922	Shield & Collar for Pedestal	
777-060	SAAR7060	Quartz Combustion Tube, SC32/132	
777-064	SAAR928	Pyrex Delivery Tube for AR7060	
777-066	SAAR7066	Quartz Liner for AR7060	
777-073	SAAR7073	Thermocouple	
778-114	SAAR8114	"O" Ring 111	1 pc.
778-116	SAAR8116	"O" Ring 256, Ballast, FP228	1 pc.
778-120	SAAR781	Neoprene Gasket, Ballast, FP228	
778-379	SAAR8380	Combustion Tube Insulator, FP528	
778-596	SAAR596	SST Filter	
778-674	SAAR674	Graphite Crucible	2g
778-891	SA778891300E	Crucible, Quartz, Large, EA	1 pc.
778-892	SA778892300E	Cover, Quartz, Large, EA	1 pc.
780-044	SAAR609	Assembly Sample Puller	
780-047	SAAR1202	Reaction Tube for DH102	
780-093	SA780093300E	Quartz Crucible, MAC400	1 pc.
780-094	SA780094300E	Quartz Crucible Cover, MAC400	1 pc.
780-379	SAAR379	Quartz Window	

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
780-831	SAAR831	"O" Ring, CHN 1000 Lance	1 pc.
780-853	SAAR853	Thermistor, Level	
780-890	SAAR7808	Large Graphite Crucible	1000 pcs.
780-892	SAAR0892	Graphite Covers for AR7808	1000 pcs.
780-899	SAAR899	Honeycomb Support, 1.25W x 1.5L FP	
780-911	SAAR910	Catalyst Heater Tube	1 pc.
780-927	SAAR927	Printer Paper Roll	
780-928	SAAR926	Printer Ribbon	
780-953	SAAR922	Pedestal CS-125	
781-098	SAAR1098	Thermistor, High Water	
781-101	SAAR101	Thermistor, Low Water	
781-165	SAAR469	SS Relay TC136 44F 6848 240 D45	
781-273	SAAR028	Ceramic Crucible Tray	
781-335	SA781335E	Quartz Boat for RC-412, RC612	1 pc.
781-371	SAAR371	Fan for CS244, 240V, 50/60Hz	
781-424	SAAR1424	"O" Ring	1 pc.
781-585	SAAR1585	Over Temp Thermocouple, RC612	
781-671	SAAR671	Ink Jet Printer Paper	
781-883	SAAR1883	Ink Jet Paper 500	
781-902	SAAR902	Ink Jet Cartridge	
782-004	SAAR928	Delivery Tube SC132/32/CR-12	
782-010	SAAR7060	Quartz Combustion Tube	
782-038	SAAR2038	Combustion Tube for RC-412	1 pc.
782-059	SAAR2059	Nickel Boat Liners, for AR1335	10 pcs.
782-079-110	SAAR9023	Ceramic Boat Stop, Hi Flow	
782-079-115	SAAR9110	Ceramic Boat Stop, Medium Flow	
782-079-120	SAAR9022	Ceramic Boat Stop, Low Flow	
782-338	SAAR2338	Moisture Removal Tube, Indicating	
782-428	SAAR2428	Quartz Catalyst Tube, RC-412	
782-534	SAAR3930	Electrode Insert Repair Service	
782-601	SAAR2601	Vacuum Cleaner	
782-703	SAAR2703	Hi Purity Outer Graphite Crucible	1000 pcs.
782-703-100	SAAR2703P	Hi Purity Outer Graphite Crucible	100 pcs.
782-719	SAAR782P	Hi Temp Graphite Crucible	200 pcs.
782-720	SAAR782	Graphite Crucible TC436	1000 pcs.
782-720	SAAR7827	Graphite Crucible, High Temp	1000 pcs.
782-720	SAAR782HD	Graphite Crucibles, Hig Temp, High Density Material	1000 pcs.
782-720S	SAAR720S	Graphite High Temp Crucible	
782-720SHD	SAAR872HD	Graphite Hi-Density Hi-Temp Crucible	
782-721	SAAR721	Hi-Temp Electrode Insert	
782-795	SAAR2795	Inner Crucible for AR2703	1000 pcs.
782-835	SAAR3930	Electrode Insert Repair Service	
782-877	SAAR877	Trim Step Brush for 782-918	
782-880	SAAR2880	Vacuum Filter for CS444	
782-918	SAAR877H	Handle Only for AR877	
782-938	SAAR938	Composite Brush for 782-942	
782-942	SAAR938H	Handle Only for AR938	
783-177	SAAR962	Heating Element for SC32/132, CR12	
783-189	SAAR3189	"O" Ring, CHN1000	1 pc.
783-190	SAAR3897	"O" Ring, CHN1000 Replaced by 783-897	1 pc.
783-193	SAAR8343	CHN1000 Furnace Element with Insulation	
783-274	SAAR278	Lance Only without Metal Ring	
783-278	SAAR278	Furnace Lance, CHN1000	
783-299	SAAR2038	Combustion Tube for RC-412	
783-301	SAAR3301	Brush End, .38 Diameter	
783-302	SAAR3302	Brush End Modified	

Leco Part No.	Säntis Part No.	Description	Qty./ pack.
783-515	SAAR515	"O" Ring	1 pc.
783-897	SAAR3897	"O" Ring CH1000	1 pc.
789-094	SAAR618	Trap Assembly, Mol Sieve, WR-12/112	
789-114	SAAR6317	Stainless Purge Tube for SC432	
789-115	SAAR789	Lance Tube Assembly, CS225, w/o Autocleaner 4.28"	
789-116	SAAR4174	Purge Tube, SC32/132	
789-118	SAAR9118	Lance Tube Assembly, 5.12 LG, CS225	
789-122	SAAR2601	Vacuum Cleaner, Specify 110/220 VAC	
789-412	SAAR9412	Catalyst Tube, TC500, RC612	
874-000-180	SA874000180W	Small Boat Loop	1 pc.
874-000-181	SA874000181W	Nickel Boat, Small for LECO AMA 245Hg	1 pc.
874-000-230	SA874000230W	Stainless Steel Boat Beam	2 pcs.
874-000-231	SA874000231W	Stainless Steel Spoiler	1 pc.
874-001-180	SA874001180W	Large Boat Loop	1 pc.
874-001-181	SA874001181W	Nickel Boat, Large for AMA 254 Hg	1 pc.

Product range: We offer thousands of products for elemental analyzers - if the product you need is not listed here contact us, as we almost certainly will have it available.

The use of Leco™ part numbers is intended for convenience only, and does not imply that the products are Leco™ origin. All Säntis Analytical products are guaranteed to be of high quality and suitable for use in the stated instrument. Leco™ is a registered trade mark.



Elementar

Productos para Analizadores Elementales

Listado de Artículos



Elementar Part No.	Säntis Part No.	Description	Qty./Pack
	SA990787	Quartz wool Disc 15mm Diameter	25 pcs.
	SA990800	Quartz wool Disc 16mm Diameter	25 pcs.
	SA990801	Quartz wool Disc 22mm Diameter	25 pcs.
	SA990802	Quartz wool Disc 19mm Diameter	25 pcs.
03 001 118	SA03001118	Sulfanilic Acid	5g
03 001 248	SA20500500	Forceps 115mm	
03 002 252	SA03002252	O-Ring 18 x 2, Red	
03 002 253/4	SA030022534	Quartz Tube	
03 002 262	SAAR049	Pt-Calalyst	50g
03 002 267	SA03002267	O-Ring 24 x 2, Green	
03 002 268	SA03002268	O-Ring 32 x 3, Red	
03 002 395	SAAR02395	O-Ring 7 x 1.5, Red	
03 002 399	SA03002399	Filter round 18x20mm for Elementar	
03 002 607	SA03002607	Shaft Sealing	
03 002 849	SAAR01040	Post Combustion Tube, (I)	
03 651 011	SA03651011	PTFE Hose 5 x 1 mm, 0.5 m	
03 653 103	SA03653103	O-Ring 9 x 2.5, Red	
03 653 118	SA03653118	O-Ring 6 x 2, Red	
03 654 505	SA03654505	O Ring Silicone Rubber	
03 654 603	SA03654603	O-Ring 11.3 x 2.4, Red	
03 654 615	SA03654615	O-Ring 8 x 2, Black	
03 654 627	SA03654627	O-Ring 20 x 2, Green	
03 654 629	SA03654629	O-Ring 25 x 2, Red	
03 654 701	SA03654701	O-Ring 4 x 4, Blue Or Red	
03 654 772	SA03654772	O-Ring 15 x 5, Green	
03 654 814	SA03654814	O-Ring 14 x 3.5, Red	
03 667 990	SA03667990	PTFE Gasket 5 mm	
03 667 991	SA03667991	Sealing Ring	
03 674 012	SA03674012	Aluminium Weighing Boats 12 x 4 x 4mm	1000 pcs.
03 674 016	SA03674016	Aluminium Weighing Boats 12 x 6 x 6mm	500 pcs.
03 679 900	SAAR168	Sansorb 20 - 40 mesh	250g
03 679 901	SAAR171	Magnesium Perchlorate	454g
03 679 905	SA99060105	Copper-Oxide, Wire	250g
03 679 907	SA990715	Quartz Chips	100g
03 679 908	SA990716	Quartz Wool	5g
03 679 910	SA03679910	High Purity Carbon Granular	80g
03 679 914	SA03679914	Lead Chromate	50g
03 679 915	SA03679915	Ceriumdioxide	100g
03 679 941	SA03679941	Square Gasket, Red	
03 679 942	SA990744B	Acetanilide	2g
03 950 005	SA20500600	Micro Spatula	
03 951 619	SA84992102	Tin Capsules	250 pcs.
03 951 620	SA84990602	Tin Capsules	250 pcs.
03 960 157	SA33840035B	Silica Gel	1kg
03 960 350	SA03960350	Nicotinic Acid Amide	1g
NEW: aquatak® Replacement for Sicapent			
	SA100050	aquatak®	50g
	SA100100	aquatak®	100g
	SA100500	aquatak®	500g
03 960 352	SA03960352	Sicapent With Indicator	500g
03 960 352	SA03960352	Sicapent With Indicator	500g
03 968 006	SAAR241	Vacuum Grease	150g
05 000 043	SA05000043	O Ring Silicone Rubber	

Elementar Part No.	Santis Part No.	Description	Qty./Pack
05 000 086	SA05000086	Al2 O3-Wool	20g
05 000 087	SA05000087	Ceramic Ash Finger	
05 000 088	SA05000088	O2-Lance, 140mm	
05 000 092	SA05000092	O-Ring 22 x 2, Black	
05 000 094	SA05000094	Gloves, Temp. Resistant	
05 000 095	SA05000095	O-Ring 22 x 2, Green	
05 000 096	SA05000096	O-Ring 18 x 2, Green	
05 000 140	SA05000140	O-Ring 9.5 x 1.8, Black	
05 000 170	SA05000170	Tin capsules for liquids 10 x 22 mm	100 pcs.
05 000 247	SA05000247	Lance-Sealing	
05 000 248	SA05000248	O-Ring 4.5 x 1.5, Red	
05 000 249	SA05000249	O-Ring 4.5 x 1.5, Green	
05 000 311	SA84994401	Tin Capsules	100 pcs.
05 000 312	SA84993001	Tin Capsules	100 pcs.
05 000 316	SAAR05316	Nitrogen/Sulfur Free Paper	500 sheets
05 000 341	SA05000341	O-Ring 6.75 x 1.78, Green	
05 000 353	SA05000353	Square Gasket, Black	
05 000 354	SA05000354	O2 - Lance, 150 mm	
05 000 359	SA05000359	O-Ring 9.19 x 2.62, Green	
05 000 368	SA05000368	O-Ring 9.25 x 1.78, Green	
05 000 369	SA05000369	O-Ring 4 x 4, Green	
05 000 371	SA05000371	O-Ring 24 x 3, Green	
05 000 372	SA05000372	O-Ring 38 x 3.5, Green	
05 000 384	SA05000384	Square gasket, Black	
05 000 385	SA05000385	O-Ring 7 x 1.5, Green	
05 000 400	SA990747B	Benzoic Acid	2g
05 000 416	SA05000416	O-Ring 7.5 x 2.5, Green	
05 000 425	SA05000425	O-Ring 5 x 1.5, Green	
05 000 429	SA76152301	Tin Foil Cups	100 pcs.
05 000 469	SA05000469	Cerium (IV) Oxide, Granular	100g
05 000 481	SA84992702	Tin Capsules	250 pcs.
05 000 658	SA05000658	O2-Lance, 170 mm, Vario EL	
05 000 699	SA99060302B	Copper, Wire	100g
05 000 814	SA05000814	Aluminium Capsules	250 pcs.
05 000 815	SA05000815	Aluminium Capsules	250 pcs.
05 000 816	SA84991506	Silver Capsules	250 pcs.
05 000 817	SA84992106	Silver Capsules	250 pcs.
05 000 880	SA05000880	Manganese (IV) Oxide	10g
05 000 900	SA990715B	Quartz Chips	50g
05 001 013	SAAR338	Gel Capsules	400 pcs.
05 001 039	SA99060305	Copper Oxide, Wires	250g
05 001 053	SA05001053	O Ring for Spring Flange	
NEW: aquatak® Replacement for Magnesium Perchlorate			
	SA100050	aquatak®	50g
	SA100100	aquatak®	100g
	SA100500	aquatak®	500g
05 001 059	SA990712B	Magnesium Perchlorate Self Indicating Granular 0.7 to 1.2mm	100g
05 001 059	SA990712B	Magnesium Perchlorate Self Indicating Granular 0.7 to 1.2mm	100g
05 001 067	SAAR01067	Glutamic Acid, Purity > 99%	250g
05 001 070	SA990700	Cobaltic Oxide, Granulated	25g
05 001 169	SA05001169	Phenylalanine	25g
05 001 540	SA76981106	Silver capsules for solids 5 x 9 mm	250 pcs.
05 001 715	SA05001715	Cotton Wool	50g

Elementar Part No.	Säntis Part No.	Description	Qty./Pack
05 001 727	SA84991502	Tin capsules for liquids 3.5 x 5.5 mm	250 pcs.
05 001 732	SA84990306	Silver capsules for liquids 2 x 5 mm	250 pcs.
05 001 928	SA84991206	Silver capsules for liquids 3.5 x 4 mm	250 pcs.
11 00 0053	SA11000053	Valve seat ring	
11 00 1150/4	SA110011504	Combustion/Reduction-Tube, Quartz	
11 00 1170/4	SA990807	Ash Finger, Quartz	
11 00 1171/4	SA110011714	Support Tube, 133 mm	
11 00 1172/4	SA110011724	Support Tube, Quartz, 65 mm	
11 00 1308/4	SA110013084	Reduction Tube, Vario EL	
11 00 1317/4	SA110013174	Ashfinger, silica 105mm	
11 00 1331/4	SA110013314	Quartz Bridge, 65 mm	
11 00 1371	SA11001371	Silver Weighing Boats 12 x 4 x 4mm	500 pcs.
11 02 0008	SA11020008	Tungsten(VI)-Oxide	60g
11 02 0017	SA11020017	Tungsten(VI)-Oxide Powder	25g
11 02 0024	SA11020024	Tungsten (VI) Oxide Powder	90g
11 02 1036/4	SA110210364	Ash Crucibles, Quartz, 60 mm	
11 02 1046	SA11021046	Silver Weighing Boats 4 x 4 x 5mm	500 pcs.
11 45 1012/4	SA114510124	Sauerstoff-Lance 190mm	
12 00 0040	SAAR02530	Tungsten, Granulated	470g
12 00 1004/4	SAAR010044	Mire Mesh, Ronde Dia. 22mm	
12 00 1005/4	SA120010054	Wire Mesh, Ronde Dia. 18 mm	
12 00 1006/4	SAAR01020	Combustion Tube RapidN	
12 00 1008/3	SAAR929	Crucible Tongs	
12 00 1009	SA12001009	Tin Boats	120 pcs.
12 00 1087/4	SA120010874	Ash Crucible, Steel	
12 01 0032	SA999990	Tin Foil, 35 x 35 mm	400 pcs.
12 01 0033	SA999991	Tin Foil, 50 x 50 mm	450 pcs.
12 01 0045	SA12010045	Tungsten Kit, Compl., Premium	
12 01 1008/4	SA12011008	Spacer	25 pcs.
12 01 1013/4	SAAR02510	Quartz Tube (Tungsten)	
12 01 1032/4	SAAR110324	Distance Ring	
12 01 1036/4	SA120110364	Combustion-, Post Combustion-Tube	
15 00 0062	SA15000062	Sulphanilamide OAS	10g
15 00 1189	SA15001189	Ash Finger - High Temperature Pyrocube	
15 00 1319	SA15001319	Ceramic Crucible	
16 00 1012	SA16001012	Drying tube 348 mm GL25	
16 00 5013	SA16005013	Consumbale Kit for approy 4000 analysis, Rapid N Cube	
22 127 245/3	SA221272453	Reduction Tube	
22 127 313/4	SA221273134	Z 623 Metal Mesh Insert	
22 131 365	SA22131365	Silver Wool	50g
22 133 210/3	SA221332103	Crack Tube GZ 384 For O2 Anal.	
22 133 211/4	SA221332114	Quartz Tube Insert	
22 133 212	SA22133212	Silver Boats	100 pcs.
22 133 213	SA22133213	Silver Boats	70 pcs.
22 133 214	SA22133214	Silver Boats	42 pcs.
22 137 418	SA22137418	Tin Boats	1000 pcs.
22 137 419	SA22137419	Tin Boats	500 pcs.
22 137 420	SA22137420	Tin Boats	250 pcs.
22 156 202/3	SA221562023	Combustion Tube	
22 160 200/4	SA221602004	Drying Tube, Compl.	
22 160 301/4	SAAR01032	U-Tube for Vario, rapid N without end caps	
22 160 301/4	SAAR01032B	U-Tube for Vario, rapid N with end caps	
23 00 1020	SA23001020	Glassy Carbon Tube Pyrocube	
23 00 1021	SA23001021	Ceramic Tube Pyrocube	
23 00 1149	SA23001149	Reaction Tube Pyrocube	
23 00 1150	SA23001150	CombustionTube Pyrocube	
25 00 0017	SAAR01007	Quartz Tube, 29 mm	

Elementar Part No.	Säntis Part No.	Description	Qty./Pack
25 00 0029	SAAR00029	Gripper Arm For CN/N	
25 00 0074	SA999992	Silver foil, 40 x 40 mm	100 pcs.
25 00 1028	SA25001028	Combustion Tube, Quartz	
25 00 1030	SAAR01030	Combustion Tube, steel	
25 00 1048	SAAR01048	Post Combustion Tube, Steel	
25 00 5042	SA25005042	Consumbale Kit - Vario MAX CHN 2000 samples Kit	
25 20 1003	SA25201003	Connecting Tube Transparent Silica 130mm (includes O-Rings)	
34 00 1339/4	SAAR22164	U-Tube, GL 18	
34 00 1069/4	SA340010694	Combustion Tube High TOC II	
35 00 2125	SA35002125	Reactor Tube liquiTOC	
50 008 449/4	SAAR084494	Wire Mesh Ronde	
50 008 467	SA50008467	Corundum Balls	40g
88 895 805	SA88895805	Aspartic Acid	100g
88 895 805	SA88895805	Aspartic Acid	



Thermo

Productos para Analizadores Elementales

Listado de Artículos



We produce capsules in sizes and materials other than those listed below. Please check our separate article list for capsules

We produce copper in sizes other than those listed below. Please check our separate article list for copper

We produce other certified standards than those listed below. Please check our separate article list for certified analytical standards

Säntis No.	Part No.	Description	Qty./pack.
Pressed Tin and Silver Capsules for solid samples			
SA76162401	252 080 00	Pressed Tin Capsules 10 x 10 mm	100
SA76980501	240 053 00	Pressed Tin Capsules 3,3 x 5 mm	100
SA76980502	240 053 00	Pressed Tin Capsules 3,3 x 5 mm	250
SA76981101	240 064 00	Pressed Tin Capsules 5 x 9 mm	100
SA76981102	240 064 00	Pressed Tin Capsules 5 x 9 mm	250
SA76981103	240 064 00	Pressed Tin Capsules 5 x 9 mm	400
SA76980705	240 054 00	Pressed Silver Capsules 4 x 6 mm	100
SA76980706	240 054 00	Pressed Silver Capsules 4 x 6 mm	250
SA76981105	240 054 05	Pressed Silver Capsules 5 x 9 mm	100
SA76981106	240 054 05	Pressed Silver Capsules 5 x 9 mm	250

Tin discs			
SA999902	252 080 18	Tin discs 35 mm, thickness 0.010 mm	100
SA999918	252 080 15	Tin discs 30 mm, thickness 0.010 mm	100
SA999976	252 080 19	Tin discs 62 mm, thickness 0.010 mm	100

Flat Bottom Smooth Tin and Silver Capsules for liquid samples			
SA84990000	240 070 00	Smooth Tin Capsules 1,06 x 5,50 mm	50
SA84990002	240 070 00	Smooth Tin Capsules 1,06 x 5,50 mm	250
SA84990301	240 088 20	Smooth Tin Capsules 2 x 5 mm	100
SA84990302	240 088 20	Smooth Tin Capsules 2 x 5 mm	250
SA84990601	240 088 10	Smooth Tin Capsules 3 x 6 mm	100
SA84990602	240 088 10	Smooth Tin Capsules 3 x 6 mm	250
SA84992101	240 072 00	Smooth Tin Capsules 3,5 x 9 mm	100
SA84992102	240 072 00	Smooth Tin Capsules 3,5 x 9 mm	250
SA84990305	240 051 00	Smooth Silver Capsules 2 x 5 mm	50
SA84990306	240 051 00	Smooth Silver Capsules 2 x 5 mm	250

Copper Reagents			
SA33835311	338 353 11	Electrolytic Copper	90 g
SA33835311B	338 353 13	Electrolytic Copper	70 g
SA33835305	338 353 05	High Quality Copper	220 g
SA33835312	338 353 12	High Quality Copper	50 g
SA33835335	338 353 12	High Quality Copper	250 g
SA99060100	338 217 10	Copper Oxide 0,7 mm	1 kg
SA99060101	338 217 10	Copper Oxide 0,7 mm	100 g
SA99060101B	338 217 10	Copper Oxide 0,7 mm	50 g
SA99060102	338 353 00	Reduced Copper Wires 0,7 mm	100 g
SA99060103	338 353 00	Reduced Copper Wires 0,7 mm	250 g
SA99060204	338 353 10	Reduced Copper Wires 0,7 mm, ultra pure	40 g

Säntis No.	Part No.	Description	Qty./ pack.
Chemical Reagents			
SA33800115	338 001 15	Glucose Anhydride	100 g
SA33835400	338 354 00	Koerbl Catalyst absorbent	5 g
SA33837530	338 375 30	Chromosorb WAW	25 g
SA33839055	338 390 55	Platinum supported on Alumina	100 g
SA33840030	338 400 30	Fluorine Adsorber	25 g
SA33840035	338 400 35	Silica Gel plus Indicator	100 g
SA990610	338 229 00	Chromium Oxide, granular	25 g
SA990610B	338 229 00	Chromium Oxide, granular	50 g
SA990620	338 354 20	Tungstic Anhydride on Alumina, granular	25 g
SA990630	338 245 00	Silvered Cobaltous Oxide, granular	25 g
SA990700	338 228 00	Cobaltous Oxide, granular	25 g
SA990701	338 238 00	Nickelized Carbon Catalyst	5 g
SA990702	338 246 00	Tungstic Anhydride, granular	25 g
SA990708	338 400 00	Oxidation Catalyst for protein analysis	40 g
SA990708C	338 400 00	Oxidation Catalyst for protein analysis	60 g
SA990709	338 375 10	Vanadium Pentoxide	1 g
SA990710	338 220 00	Carbon Dioxide Absorbent	100 g
SA990711	338 375 40	Lead Oxide, red	1 g
NEW: aquatak® Replacement for Magnesium Perchlorate			
SA100050		aquatak®	50g
SA100100		aquatak®	100g
SA100500		aquatak®	500g
SA990712	338 219 00	Magnesium Perchlorate, granular	100 g
SA990712B	338 219 00	Magnesium Perchlorate, granular, Self Indicating	100 g
SA990713	338 018 01	Molecular Sieve 3A, pellets, gas purification	100 g
SA990714	338 250 00	Nickel Wool	2 g
SA990715	338 223 00	Quartz Chips, granular	100 g
SA990715B	338 223 00	Quartz Chips, granular	50 g
SA990716	338 222 00	Quartz Wool	5 g
SA990716A	338 222 05	Quartz Wool	15 g
SA990716B	338 222 00	Quartz Wool	50 g
SA990717	-	Soda Lime, 1.7 - 5.0 mm	100 g
SA990717A	-	Soda Lime, 1.0 - 2.5 mm	100 g
SA990717B	-	Soda Lime, 1.7 - 5.0 mm	200 g
SA990717C	-	Soda Lime, 1.7 - 5.0 mm, Self Indicating	100 g
SA990717D	-	Soda Lime, 1.0 - 2.5 mm, Self Indicating	100 g
SA999985	338 220 00	Ascarite	100 g

Seperation Columns			
SA26007820	260 078 20	Packed Column Nitrogen, 0,5 mt PTFE	1
SA26007920	260 079 20	Multiseperation Column SST, 2mt, 6x5	1
SA26008207	260 082 07	Packed Column CHN 2 mt PTFE	1
SA26008208	260 082 08	Packed Column 2 mt PTFE	1
SA26008215	260 082 15	Multisep.Column CHNS 2mt 6x4, PTFE	1

Säntis No.	Part No.	Description	Qty./ pack.
SA26008220	260 082 20	Multiseparation Column PTFE 2mt 6x5mm	1
SA26008245	260 082 45	Packed Column CHN/4 mt, SST	1
SA26070004	260 700 04	Teflon Column for Protein Determination	1
SA990720	260 082 06	Packed Column N2/CO2/H2O/SO2 2 mt	1
SA990721	260 082 00	Packed Column N2/CO2/H2O 2 mt	1
SA990722	260 082 05	Packed Column N2/CO2 3 mt	1
SA990723	260 078 00	Packed Column SO2 0,8 mt	1
SA990724	260 079 00	Packed Column O2 1 mt	1

Certified Organic Analytical Standards			
SA33800110	338 001 10	EDTA Standard	100 g
SA33835430	338 354 30	Imidazole Standard	2 g
SA33840001	338 400 01	Urea Standard	2 g
SA33840002	338 400 02	Urea Standard	10 g
SA33840018	338 400 18	L-Cystine Standard	1 g
SA33840018B	338 400 18	L-Cystine Standard	2 g
SA33840019	338 400 19	Nicotinamide Standard	2 g
SA33840022	338 400 22	Aspartic Acid	2 g
SA33840023	338 400 23	Aspartic Acid Standard	10 g
SA33840024	338 400 24	Pasta Reference Material	10 g
SA33840027	338 400 27	Pasta Reference Material for S Determination	2g
SA33840028	338 400 28	Pasta Reference Material	100 g
SA33840025	338 400 25	Soil Reference Material	5 g
SA33840026	338 400 26	Soil Reference Material CNS	5 g
SA990744	338 367 00	Acetanilide Standard	1 g
SA990744B	338 367 00	Acetanilide Standard	2 g
SA990746	338 244 00	Atropine Standard	1 g
SA990746B	338 244 00	Atropine Standard	2 g
SA990747	338 225 00	Benzoic Acid Standard	1 g
SA990747B	338 225 00	Benzoic Acid Standard	2 g
SA990748	338 224 00	Cyclohexanone 2 4 Dinitr. Hydr.	1 g
SA990750	338 251 00	Sulphanilimide Standard	1 g
SA990750B	338 251 00	Sulphanilimide Standard	2 g
SA990751	338 352 20	Methionine	10 g
SA990751B	338 352 20	Methionine	25 g
SA990752	338 352 10	BBOT Standard	1 g
SA990752B	338 352 10	BBOT Standard	2 g
SA999954	338 352 40	Lubricant Standard	2 g

Quartz- and Metallic Tubes			
SA46820070	468 200 70	Empty Quartz Reactor Tube, 18 mm (Flash EA 1112)	2
SA990729	468 200 60	Empty Quartz Reactor Tube, 18 mm (NA1500)	2
SA990731	-	Opaque quartz tube 18 x 450 mm (NA1500, EA1108)	2
SA990735	468 016 00	Empty Quartz Reactor Tube, 20 mm (NA 1500)	2
SA46820010	468 200 50	Empty SST Combustino tube 18mm (NA2000)	1
SA46820075	468 200 75	Empty Inconel Reactor Tube, 25mm OD (HPAR)	1
SA990725	468 009 00	Empty Quartz Combustion Tube (1106)	5

Säntis No.	Part No.	Description	Qty./ pack.
SA990726	468 010 00	Empty Quartz Reduction Tube (1106)	2
SA25204510	252 045 10	Quartz Insert, 210 mm, 13 mm	1
SA999980	252 045 10	Quartz Insert, 235 mm, 13 mm	1
SA999938	-	Quartz Insert 210 mm, 12mm (NA1500/1108/1110)	5
SA25205001	252 050 10	Inconel Crucible, 13mm with slots	1
SA25205011	252 050 11	Inconel crucible 18mm OD (Flash EA 1112) regular slits	1
SA25205011B	252 050 11	Inconel crucible 18mm OD (Flash EA 1112) less slits	1
SA27204505	272 045 05	Ceramic Disk	20
SA990740	281 131 00	Water Trap W/OR small caps	1
SA990740B	-	Plug and septa for trap tube for H2O	1
SA990740C	231 010 08	Small trap tube for H2O without plug and septa	1
SA990741	281 130 90	Great trap tube for CO2 with plug and septa	1
SA990741C		Great trap tube for CO2 without plug and septa	1
SA990733	468 020 05	Prepacked Quartz Reactor Tube (EA1108, NA1500) Electrolytic Copper/Tungstic Anhydride, CHNS	1
SA46802007	468 020 07	Prepacked Quartz Reactor Tube(Flash EA1112) Electrolytic Copper/Copper Oxide	1
SA46802008	468 020 08	Prepacked HPAR Reactor Tube (Flash EA1112) with oxidation catalyst	1
SA999977	468 020 13	Prepacked Quartz Reactor Tube (Flash EA 1112) Chromium-/Silvered cobaltous oxide/high purity copper	1
SA46820051	468 200 51	Prepacked SS Reactor Tube (NA1500/2, EA1108) Electrolytic Copper/Tungstic Anhydride	1
SA999967	-	Prepacked Quartz Reactor Tube (Flash EA 1112) Chromium oxide/Silvered Cobaltous oxide	1

O-Rings and Seals			
SA29000051	290 000 51	O-Ring lubricant	1
SA29003603	290 036 03	O-ring Silicon for H2O-Trap	1
SA29003635	290 036 35	Sealing o-ring for Pyrex tube	2
SA29010386	290 103 86	Viton O-Ring for upper and button connection	5
SA29010684	290 106 84	Silicon O-Ring for HPAR 4000 crucible	5
SA29020640	290 206 40	Bottom O-Ring for 18mm OD	5
SA29020649	290 206 49	Bottom O-Ring for HPAR Reactor	5
SA29020682	290 206 82	Top O-Ring for HPAR Reactor	5
SA29022910	290 229 10	Top O-Ring for 18mm OD	5
SA29030304	290 303 04	o-ring viton 1/8"	10
SA29032940	290 329 40	O-Ring Cone for AS128 auto sampler	1
SA34301512	343 015 12	Shaft for autosampler MAS200R	1
SA34301512R		O-Ring, green, for Shaft	3
SA34301512WZ		Tool to replace the O-Ring (SA34301512R)	1
SA990727	290 303 31	O-Ring Viton Top	5
SA990728	290 305 67	Bottom O-Ring	10
SA990732	290 329 10	O-Ring Viton	10
SA990742	290 136 03	Seals for Traps	2

Gas Fittings and Tubing			
SA29013486	290 134 86	Ferrule graphite, i.d. 0.8	2
SA29013487	290 134 87	Ferrule graphite, i.d. 0.45	2

Säntis No.	Part No.	Description	Qty./ pack.
SA29013488	290 134 88	Ferrule graphite, i.d. 0.35	2
SA29034044	290 340 44	Tubing Olive 1,6 mm, aluminium	10
SA29034045	290 340 45	Tubing Olive 2 mm, aluminium	10
SA29034151	290 341 51	Tubing Olive 4 mm, aluminium	10
SA29034245	290 342 45	Tubing Olive 2 mm, stainless steel	10
SA29034545	290 345 45	Tubing Olive 2 mm, PTFE	10
SA34714100	347 141 00	Connector union 2mm: Set with Olive, Nut and Union	1
SA34715401	347 154 01	Column Fitting, 6 - 2mm	1
SA34734100	347 341 00	Connector union 6MB-6MB	5
SA34736100	347 361 00	Tubing Connector Block, 4 way, 2mm	5
SA34738103	347 381 03	T-Connection Brass, 3x6 MB	5
SA35010100	350 101 00	Column Fitting, 18 - 2mm	1
SA35034108	350 341 08	Rubber Holder SST 6MB	5
SA35040403	350 404 03	Stainless steel nuts 6MB	10
SA35045100	350 451 00	Blanking nut 6MB	10
SA42001501	420 015 01	PTFE tube 2x1mm	5 mt
SA42001601	420 016 01	PTFE Tubing 4x2mm	5 mt
SA42001701	420 017 01	PTFE Tubing 6x4mm	5 mt
SA42010000	420 100 00	Tube PE- 5x3mm	5 mt
SA39102500	391 025 00	Stainless steel tube 2mm	5 mt
SA590200		Cutter for SST-tube	1
SA999922		Ferrule graphite, 1/4" to 1/16"	10

Miscellaneous and Spare Parts			
SA19004157	190 041 57	Cone for AS128 auto sampler	1
SA20500310	205 003 10	Sealing device for solid samples (tweezer pencil+plate)	1
SA999960B	205 003 10	Tweezer pencil	1
SA999960	205 003 10	Plate with 4 holes	1
SA20500500	205 005 00	Forceps	1
SA20500600	205 006 00	Spatula for container filling (liquids and solids)	1
SA20503002	205 030 02	Sealing device for liquid samples	1
SA24806400	248 064 00	Thermal paper for Printer	10
SA31302100	313 021 00	Membrane for bi-pneumatique valve	10
SA34301512	343 015 12	Shaft for autosampler MAS200R	1
SA35460015B	354 500 10	Heating Element NA1500/1108	1
SA35460021	354 600 21	Heating Element for Flash EA1112	1
SA25117400	251 174 00	MAS200R Autosampler	1
SA37001101	370 011 01	Slider for AS 200 LS auto sampler	1
SA40311095	403 110 95	Temperature Sensor Pt100 for Oven	1
SA40320000	403 200 00	Thermocouple	1
SA40320005	403 200 05	Thermocouple Pt/Pt-Rh for Flash EA1112	1
SA40320006	403 200 06	Thermocouple for Flash HT	1
SA40510100	405 101 00	Micro Valve 4way	1
SA40600300	406 003 00	Electric fan	1
SA999931	334 062 00	Nickel Sheet 130 x 130 mm	1

TC/EA Consumables

Säntis No.	Part No.	Description	Qty./ pack.
SA1041480	1041480	Ceramic Tube, Length 320mm, ID 0.5mm, OD 1.5mm	1
SA1117332	1117332	Graphite Crucibles	1
SA1117400	1117400	Glassy Carbon Granulate	50 g
SA1121030	1121030	GC Separation Column for TC/EA, 5A, 0,6 m Mesh 80-100, ID 4mm, OD 1/4"	1
SA1121120	1121120	O-ring for TC/EA	1
SA1121310	1121310	Glassy Carbon Reactor Tube	1
SA1121310A		Glassy Carbon Reactor Tube 450mm	1
SA1121310B		Glassy Carbon Reactor Tube 356mm	1
SA1121320	1121320	Ceramic Tube, OD 17mm, ID 13mm, 470mm	1
SA1121320MOD		Ceramic Tube, OD 18mm, ID 14mm, 470mm	1
SA1121330	1121330	Bottom reactor connection screw for TC/EA	1
SA1121470	1121470	Ceramic tube with thermo element for TC/EA	1
SA1156890	1156890	Ceramic Tube, OD 17mm, ID 13mm, 450mm	1
SA990786	1117430	Silver Wool	10 g
SA990786B	1117430	Silver Wool	2 g
SA990786C	1117430	Silver Wool	50 g

Product range: We offer thousands of products for elemental analyzers - if the product you need is not listed here

The use of Thermo™ part numbers is intended for convenience only, and does not imply that the products are Thermo™ origin. All Säntis



Eurovector

Productos para Analizadores Elementales

Listado de Artículos



Article List for Elemental Analyzer Eurovector™

We produce capsules in sizes and materials other than those listed below. Please check our separate article list for capsules

We produce copper in sizes other than those listed below. Please check our separate article list for copper

We produce other certified standards than those listed below. Please check our separate article list for certified analytical standards

Säntis No.	Part No.	Description	Qty./ pack.
Pressed Tin and Silver Capsules for solid samples			
SA76162401	E12010	Pressed Tin Capsules 10 x 10 mm	100
SA76980501	E12003	Pressed Tin Capsules 5,0 x 3,3 mm	100
SA76980502	E12001	Pressed Tin Capsules 5,0 x 3,3 mm	250
SA76981101	E12005	Pressed Tin Capsules 5 x 9 mm	100
SA76981102	E12007	Pressed Tin Capsules 5 x 9 mm	250
SA76980505	E12050	Pressed Silver Capsules 5,0 x 3,3 mm	100
SA76980506	E12050	Pressed Silver Capsules 5,0 x 3,3 mm	250
SA76981105	E12055	Pressed Silver Capsules 5 x 9 mm	100
SA76981106	E12055	Pressed Silver Capsules 5 x 9 mm	250

Flat Bottom Smooth Tin and Silver Capsules for liquid samples			
SA84990300	E12020	Smooth Tin Capsules 2 x 5 mm	50
SA84990301	E12020	Smooth Tin Capsules 2 x 5 mm	100
SA84990302	E12020	Smooth Tin Capsules 2 x 5 mm	250
SA84990601	E12025	Smooth Tin Capsules 2.9 x 6 mm	100
SA84990602	E12025	Smooth Tin Capsules 2.9 x 6 mm	250
SA84992101	E12030	Smooth Tin Capsules 3,5 x 9 mm	100
SA84992102	E12030	Smooth Tin Capsules 3,5 x 9 mm	250

Copper Reagents			
SA99060100	E10108	Copper Oxide 0,7 mm	1 kg
SA99060101	E10108	Copper Oxide 0,7 mm	100 g
SA99060101B	E10108	Copper Oxide 0,7 mm	50 g
SA99060102	E10121	Reduced Copper Wires 0,7 mm	100 g
SA99060103	E10122	Reduced Copper Wires 0,7 mm	250 g
SA99060204	E10130	Reduced Copper Wires 0,7 mm, ultra pure	40 g

Chemical Reagents			
SA33837530	E10504	Chromosorb WAW	25 g
SA990610	E10101	Chromium Oxide, granular	25 g
SA990610B	-	Chromium Oxide, granular	50 g
SA990620	E10156	Tungstic Anhydride on Alumina, granular	25 g
SA990630	E10160	Silvered Cobaltous Oxide, granular	25 g
SA990700	E10105	Cobaltous Oxide, granular	25 g
SA990701	E10146	Nickelized Carbon Catalyst	5 g
SA990702	E10152	Tungstic Anhydride, granular	25 g
SA990708	E10148	Oxidation Catalyst for protein analysis	40 g
SA990708C	-	Oxidation Catalyst for protein analysis	60 g
SA990709	E10184	Vanadium Pentoxide	1 g
SA990710	E10503	Carbon Dioxide Absorbent	100 g

Säntis No.	Part No.	Description	Qty./ pack.
SA990711	E10180	Lead Oxide, red	1 g
NEW: aquatak® Replacement for Magnesium Perchlorate			
SA100050		aquatak®	50g
SA100100		aquatak®	100g
SA100500		aquatak®	500g
SA990712	E10507	Magnesium Perchlorate, granular	100 g
SA990712B	-	Magnesium Perchlorate, granular, Self Indicating	100 g
SA990713	E10508	Molecular Sieve 3A, pellets, gas purification	100 g
SA990714	E10510	Nickel Wool	2 g
SA990715	E10513	Quartz Chips, granular	100 g
SA990715B	E10513	Quartz Chips, granular	50 g
SA990716	E10514	Quartz Wool	5 g
SA990716A	E10514	Quartz Wool	15 g
SA990716B	E10514	Quartz Wool	50 g
SA990717	-	Soda Lime, 1.7 - 5.0 mm	100 g
SA990717A	-	Soda Lime, 1.0 - 2.5 mm	100 g
SA990717B	-	Soda Lime, 1.7 - 5.0 mm	200 g
SA990717C	-	Soda Lime, 1.7 - 5.0 mm, Self Indicating	100 g
SA990717D	-	Soda Lime, 1.0 - 2.5 mm, Self Indicating	100 g
SA990784	E10504	Chromosorb W 30 - 60 mesh	10 g
SA990784B	E10504	Chromosorb W 120-140 mesh	10 g
SA990817	E10142	Nickel oxide pure	50 g
SA999931	-	Nickel Sheet 130 x 130 mm	1

Tubes			
SA999940	E12509	Quartz Empty Reactor Tube, 450mm, 18/6mm	2
SA999940B	E12510	Quartz Empty Reactor Tube, 450mm, 18/18mm	2
SA999989	E12512	Quartz Empty Reactor Tube for (H) in liquids (Modified)	1
SA46822118	HE46830020	Quartz Empty Reactor Tube, 435mm, 21/18mm	2
SA1121320	E12507	Ceramic Tube, 470mm, OD 17mm, ID 13mm	1
SA999994	HE468110	Ceramic Tube, 472mm, OD 19mm, ID 15mm	1
SA999988	E12529	Glassy Carbon Reactor Tube, 400mm, OD 12, ID 9mm	1
SA999996	E12531	Quartz Insert, 220mm, 13mm	1
SA25204510	E12531	Quartz Insert, 210mm, 13 mm	1
SA999980	-	Quartz Insert, 235mm, 13mm	1
SA999980B	-	Quartz Insert, 235mm, 13mm (one end closed)	1
SA25205010	-	Refractory metallic Crucible, 250 mm, OD 13 mm	1
SA46850020	HE46850020	Ceramic Insert, 235mm, OD 17mm, ID 13mm	2
SA990740	E13502	Small trap tube for H ₂ O, with plug and septa	1
SA990741	E13508	Great trap tube for CO ₂ , with plug and septa	1

Packed Reactors			
SA13001	E13001	Prepacked Reactor Tube 18/18mm, CHNS Tungstic Ox/Al, Red Copper, Quartz chips	1
SA13005	E13005	Prepacked Reactor Tube 18/18mm, CN/N Chrom. Ox, Silver Cobalt ox, Quartz chips	1

Säntis No.	Part No.	Description	Qty./ pack.
SA13010	E13010	Prepacked Reactor Tube 18/18mm, CN/N Quartz Chips, Red. Copper, Quartz Chips	1
SA46820699B	E13025	Prepacked Reactor Tube 18/18mm, Oxygen Pyrolysis Nickel Wool, Nickel. Carbon, Quartz Chips	1
SA13030	E13030	Prepacked Reactor Tube 18/18mm, CHN Chrom. Ox, Red. Copper, Silv. Cob. Ox, Quartz Chips	1
SA13032	E13032	Prepacked Reactor Tube 18/6mm, N-Protein	1
SA13036	E13036	Prepacked Reactor Tube Combust. 18/6mm, CHN/CN/N Chrom ox granul(E10101), Silver cobalt ox,	1
SA13037	E13037	Prepacked Reactor Tube Reduction 18/6mm, CHN/CN/N Copper ox (E10108), Red. Copper (E10120)	1
SA46830500B	E13040	Prepacked Reactor Tube Single 18/6mm, CHNS/CNS/S Tungstic Ox, Red. Copper for S	1
SA46830500	E13040	Prepacked Reactor Tube Single 18/6mm, CHNS/CNS/S Tungstic Ox, Red. Copper	1
SA13041	E13041	Prepacked Reactor Tube Single 18/6mm, CHN/CN/N Chromium Ox, Red. Copper, Silver Cob. Ox,	1
SA46820699	E13043	Prepacked Reactor Tube 18/6mm Single, Oxygen Nickel. Carbon, Quartz Chips	1
SA46830600	HE46830600	Prepacked Reactor Tube, 435mm, CHNS Copper for S / Tungstic Ox / Al	1
SA46820995	HE46820995	Prepacked Reactor Tube, 18/6 mm, CHN, CN, N WO3/Al and CoO/Ag	1
Seperation Columns			
SA990720	E11501	Packed Column N2/CO2/H2O/SO2 2 mt	1
SA990721	E11505	Packed Column N2/CO2/H2O 2 mt	1
SA26008208	E11507	Packed Column 2 mt PTFE	1
SA990722	E11510	Packed Column N2/CO2 3 mt	1
SA990723	E11515	Packed Column SO2 0,8 mt	1
SA990724	E11520	Packed Column O2 1 mt	1
SA999995	E11521	Seperation Column O/H, 1.5m, SST	1
SA26008201	HE26008201	CHNS Column 1m, PTFE	1
SA26008245		Packed Column CHN/4 mt, SST	1
SA26070004	HE26070004	Teflon Column for Protein Determination	1
SA999986	HE26070500	Column for CN/ 60 cm	1

Certified Organic Analytical Standards			
SA990744	E11001	Acetanilide Standard	1 g
SA990744B	E11001	Acetanilide Standard	2 g
SA990762	E11002	4-Fluorobenzoic Acid	1 g
SA33840001	E11003	Urea Standard	2 g
SA990745	E11004	Anthracène	1 g
SA502092	E11005	EDTA	50 g
SA990753	E11007	Atropine Sulphate	1 g
SA990746	E11008	Atropine Standard	1 g
SA990746B	E11008	Atropine Standard	2 g
SA33840018	E11009	L-Cystine Standard	1 g
SA990747	E11016	Benzoic Acid Standard	1 g
SA990747B	E11016	Benzoic Acid Standard	2 g

Säntis No.	Part No.	Description	Qty./ pack.
SA990748	E11020	Cyclohexanone 2 4 Dinitr. Hydr.	1 g
SA990749	E11024	1Cloro 2 4-Dinitrobenzen	1 g
SA990750	E11028	Sulphanilimide Standard	1 g
SA990750B	E11028	Sulphanilimide Standard	2 g
SA990751	E11026	Methionine	10 g
SA990751B	E11026	Methionine	25 g
SA990752	E11012	BBOT Standard	1 g
SA990752B	E11012	BBOT Standard	2 g

O-Rings and Seals			
SA29013520	E13520	O-Ring viton for 10/11mm OD	10
SA990728	E13525	O-Ring Viton, 6 mm	10
SA29020640	E13530	Bottom O-Ring for 18mm OD	5
SA29013535	E13535	O-Ring EVR for holder	2
SA29013536	E13536	O-Ring EVR for reactor 18mm OD	1
SA29013537	E13537	O-Ring EVR for bottom quick fit 18mm	1
SA29013538	E13538	O-Ring for Ceramic Reactor 17mm	1
SA29013539	E13539	O-Ring EVR for bottom quick fit 17mm	1
SA29013540	E13540	O-Ring EVR for bottom quick fit 6mm	1
SA29003603	HE29003603	O-ring Silicon for H2O-Trap	1
SA29003635	HE29003635	Sealing o-ring for Pyrex tube	2

Gas Fittings and Tubing			
SA29013486	290 134 86	Ferrule graphite, i.d. 0.8	2
SA29013487	290 134 87	Ferrule graphite, i.d. 0.45	2
SA29013488	290 134 88	Ferrule graphite, i.d. 0.35	2
SA29034044	290 340 44	Tubing Olive 1,6mm, aluminium	10
SA29034045	290 340 45	Tubing olive 2mm, aluminium	10
SA29034245	290 342 45	Tubing olive 2mm, stainless steel	10
SA29034545	290 345 45	Tubing olive 2mm, PTFE	10
SA34714100	347 141 00	Connector union 2mm	1
SA34715401	347 15 01	Column Fitting, 6 - 2mm	1
SA34734100	347 341 00	Connector union 6MB-6MB	5
SA34736100	347 361 00	Tubing Connector Block, 4 way, 2mm	5
SA34738103	347 381 03	T-Connection Brass, 3x6 MB	5
SA35010100	350 101 00	Column Fitting, 18 - 2mm	1
SA35034108	350 341 08	Rubber Holder SST 6MB	5
SA35040403	350 404 03	Stainless steel nuts 6MB	10
SA35045100	350 451 00	Blanking nut 6MB	10
SA42001501	420 015 01	PTFE tube 2x1mm	5 mt
SA42001601	420 016 01	PTFE Tubing 4x2mm	5 mt
SA42001701	420 017 01	PTFE Tubing 6x4mm	5 mt
SA42010000	420 100 00	Tube PE- 5x3mm	5 mt
SA999922		Ferrule graphite, 1/4" to 1/16"	10

Miscellaneous and Spare Parts			
SA20500310	205 003 10	Sealing device for solid samples (tweezer pencil+plate)	1

Säntis No.	Part No.	Description	Qty./ pack.
SA999960B	205 003 10	Tweezer pencil	1
SA999960	205 003 10	Plate with 4 holes	1
SA20500500	205 005 00	Forceps	1
SA20500600	205 006 00	Spatula for container filling (liquids and solids)	1
SA20503002	205 030 02	Sealing device for liquid samples	1
SA39102500	391 025 00	Stainless steel tube 2mm	5 mt

Isotope Grade consumables			
SA990610	E10101-IG	Chromium Oxide, granular	25 g
SA99060204	E10130-IG	Reduced Copper Wires 0,7 mm, ultra pure	40 g
SA990702	E10152-IG	Tungstic Anhydride, granular	25 g
SA999987	E10521-IG	Glassy Carbon Granulate	50 g
SA33840001	E11003-IG	Urea Standard	2 g
SA990752	E11012-IG	BBOT Standard	1 g
SA999995	E11521-IG	Seperation Column O/H, 1.5m, SST	1
SA1121320	E12507-IG	Ceramic Tube, 470mm, OD 17mm, ID 13mm	1
SA999989	E12512-IG	Modified Quartz Reactor Tube for (H) in liquids	1
SA999988	E12529-IG	Glassy Carbon Reactor Tube, 400mm, OD 12, ID 9mm	1

Product range: We offer thousands of products for elemental analyzers - if the product you need is not listed here contact us, as we almost certainly will have it available.

The use of Eurovector™ part numbers is intended for convenience only, and does not imply that the products are Eurovector™ origin. All Säntis Analytical products are guaranteed to be of high quality and suitable for use in the stated instrument. Eurovector™ is a registered trade mark.



MATERIALES DE REFERENCIA CERTIFICADOS PARA ANALIZADORES ELEMENTALES

Listado de Artículos



EQUILAB

SANTIS
analytical

Price List for Certified Analytical Standards

valid from 1st May 2010

These purified organic compounds are suitable for use as certified working standards in Elemental Analysis and are supplied with batch Certification of Analysis.

Sántis No.	Standard	C %	H %	N %	O %	S %	Qty./ pack.
SA33800110	EDTA Standard	41.09	5.52	9.58	43.80	-	100 g
SA502092	EDTA Standard	41.09	5.52	9.58	43.80	-	50 g
SA33800115	Glucose Anhydride	39.99	6.72	-	53.28	-	100 g
SA33802015	P-Nitroanilin	52.17	4.38	20.28	23.17	-	1g
SA33835430	Imidazole	52.93	5.92	41.15	-	-	2 g
SA33840001	Urea	20.00	6.71	46.64	26.64	-	2 g
SA33840002	Urea	20.00	6.71	46.64	26.64	-	10 g
SA33840018	L-Cystine	29.99	5.03	11.66	26.63	26.69	1 g
SA33840018B	L-Cystine	29.99	5.03	11.66	26.63	26.69	2 g
SA33840019	Nicotinamide	59.01	4.95	22.94	13.10	-	2 g
SA33840023	Aspartic Acid	36.09	5.30	10.52	48.08	-	10 g
SA33840024	Pasta Reference Material	-	-	2.09	-	-	10 g
SA33840025	Soil reference material	3.50	-	0.37	-	-	5 g
SA33840026	Soil reference material CNS	3.50	-	0.37	-	0.06	5 g
SA33840027	Pasta Reference Material	-	-	-	-	0.14	2 g
SA990744	Acetanilide	71.09	6.71	10.36	11.84	-	1 g
SA990744B	Acetanilide	71.09	6.71	10.36	11.84	-	2 g
SA990745	Anthracene	94.34	5.66	-	-	-	1 g
SA990746	Atropine	70.56	8.01	4.84	16.59	-	1 g
SA990746B	Atropine	70.56	8.01	4.84	16.59	-	2 g
SA990747	Benzoic Acid	68.85	4.95	-	26.20	-	1 g
SA990747B	Benzoic Acid	68.85	4.95	-	26.20	-	2 g
SA990748	Cyclohexanone 2-4 / DNPH	51.79	5.07	20.14	23.00	-	1 g
SA990749	1-Chloro-2,4-Dinitrobenzene	35.58	1.49	13.83	31.60	-	1 g
SA990750	Sulphanilamide	41.85	4.68	16.26	18.58	18.62	1 g
SA990750B	Sulphanilamide	41.85	4.68	16.26	18.58	18.62	2 g
SA990751	Methionine	40.25	7.43	9.39	21.44	21.49	10 g
SA990751B	Methionine	40.25	7.43	9.39	21.44	21.49	25 g
SA990752	BBOT	72.52	6.09	6.51	7.43	7.44	1 g
SA990752B	BBOT	72.52	6.09	6.51	7.43	7.44	2 g
SA990753	Atropine Sulphate	58.77	7.25	4.03	25.33	4.61	1 g
SA990754	4-Bromobenzoic Acid	41.82	2.51	-	15.92	-	1 g
SA990755	2-Iodobenzoic Acid	33.90	2.03	-	15.92	-	1 g
SA990756	Sulphanilic Acid	41.61	4.07	8.09	27.72	18.50	1 g
SA990772	Diphenyl OAS Standard	93.46	6.54	-	-	-	1 g
SA990775	Ammonium Sulphate standard	-	6.10	21.20	48.43	24.27	5 g
SA990775A	Ammonium Sulphate standard	-	6.10	21.20	48.43	24.27	10 g
SA990812	Phenantrene	93.55	5.62	0.18	-	0.46	1 g
SA999946	Caffeine	49.48	5.19	28.85	16.48	-	1 g
SA999954	Lubricant	72.36	11.44	0.90	-	3.06	2 g

We can supply other standards with certificates of analysis upon request. So please don't hesitate to contact us.



Leco

STANDARDS

Listado de Artículos



Article List for Leco/Eltra Standards

Säntis No.	Leco No.	Size	C	H	N	O	S	Qty./ pack.
Carbon / Sulfur Ring Standards								
SAAR869		1 g	0.0190 %	-	-	-	0.0140 %	454 g
SAAR870	501-501	1 g	0.0320 %	-	-	-	0.0330 %	454 g
SAAR871	501-502	1 g	0.0510 %	-	-	-	0.0340 %	454 g
SAAR872	501-503	1 g	0.1800 %	-	-	-	0.0190 %	454 g
SAAR873	501-504	1 g	0.4680 %	-	-	-	0.0109 %	454 g
SAAR874	501-505	1 g	0.7230 %	-	-	-	0.0220 %	454 g
SAAR875	501-506	1 g	0.8290 %	-	-	-	0.0067 %	454 g
SAAR876		1 g	0.5460 %	-	-	-	0.0082 %	454 g
SAAR878		1 g	0.2710 %	-	-	-	0.0190 %	454 g

Carbon / Sulfur Pin Standards								
SAAR511		0.5 g	0.0810 %	-	-	-	0.2400 %	250 g
SAAR512		0.5 g	0.0150 %	-	-	-	0.2620 %	250 g
SAAR881	501-674	1 g	0.0140 %	-	-	-	0.0045 %	454 g
SAAR882	501-675	1 g	0.0630 %	-	-	-	0.0230 %	454 g
SAAR883		1 g	0.1900 %	-	-	-	0.0220 %	454 g
SAAR884	501-678	1 g	0.4410 %	-	-	-	0.0110 %	454 g
SAAR885	501-679	1 g	0.8050 %	-	-	-	0.0220 %	454 g
SAAR886		1 g	0.7380 %	-	-	-	0.0240 %	454 g
SAAR888		1 g	0.0820 %	-	-	-	0.2430 %	454 g
SAAR889		1 g	1.0570 %	-	-	-	0.0215 %	454 g
SAAR890	501-676	1 g	0.0990 %	-	-	-	0.0182 %	454 g
SAAR891		1 g	0.0380 %	-	-	-	0.0330 %	454 g
SAAR892		1 g	0.0049 %	-	-	-	0.0020 %	454 g
SAAR893		1 g	1.5100 %	-	-	-	0.0130 %	454 g
SAAR894		1 g	0.5440 %	-	-	-	0.0096 %	454 g
SAAR895	501-677	1 g	0.2470 %	-	-	-	0.0210 %	454 g
SAAR896		1 g	0.1680 %	-	-	-	8.7000 %	454 g

Carbon / Sulfur / Nitrogen Chip Standards								
SAAR946			0.0029 %	-	0.0019 %	-	0.0043 %	150 g
SAAR947			0.3890 %	-	121 PPM	-	0.0940 %	150 g
SAAR948			0.5030 %	-	93 PPM	-	0.0770 %	150 g
SAAR949			0.0137 %	-	860 PPM	-	0.0270 %	150 g
SAAR950			0.0550 %	-	55 PPM	-	0.0090 %	150 g
SAAR951			0.1470 %	-	0.0116 %	-	0.0270 %	150 g
SAAR952			0.4700 %	-	104 PPM	-	0.0210 %	150 g
SAAR953			0.9270 %	-	57 PPM	-	0.0110 %	150 g
SAAR954			0.4240 %	-	84 PPM	-	0.0310 %	150 g
SAAR955			0.2040 %	-	91 PPM	-	0.0284 %	150 g
SAAR956			0.4980 %	-	0.0097 %	-	0.3010 %	150 g
SAAR957			0.1540 %	-	30 PPM	-	0.1120 %	150 g
SAAR958			0.0150 %	-	0.0680 %	-	0.0230 %	150 g

Säntis No.	Leco No.	Size	C	H	N	O	S	Qty./ pack.
SAAR960			0.1330 %	-	0.0350 %	-	0.0013 %	150 g
SAAR961			0.0220 %	-	0.0410 %	-	0.0093 %	150 g

Cast Iron thru wrought Iron Standards								
SAAR299			1.0500 %	-	-	-	0.0110 %	150 g
SAAR300			3.3200 %	-	-	-	0.0140 %	150 g
SAAR301	501-105		2.0200 %	-	-	-	0.0180 %	150 g
SAAR302			3.2900 %	-	-	-	0.0460 %	150 g
SAAR303			3.9400 %	-	-	-	0.0220 %	150 g
SAAR304	501-083		4.2300 %	-	-	-	0.0140 %	150 g
SAAR305			2.7900 %	-	-	-	0.0120 %	150 g
SAAR306			3.5000 %	-	-	-	0.0190 %	150 g
SAAR307	501-123		6.1500 %	-	-	-	0.0041 %	150 g
SAAR309			3.7000 %	-	-	-	0.0240 %	150 g
SAAR310			3.3500 %	-	-	-	0.0850 %	150 g
SAAR315			3.0900 %	-	-	-	0.1110 %	150 g
SAAR318			3.0800 %	-	-	-	0.0370 %	150 g
SAAR319			3.2800 %	-	-	-	0.0550 %	150 g
SAAR322	501-024		3.4100 %	-	-	-	0.0500 %	250 g
SAAR323			3.8400 %	-	-	-	0.0640 %	250 g
SAAR673	501-673 502-231		0.0004 %	-	-	-	0.0017 %	454 g

Oxygen / Nitrogen in Steel Pins								
SAAR644	501-644	1 g	-	-	157 PPM	3 PPM	-	100 pcs.
SAAR645		1 g	-	-	60 PPM	52 PPM	-	100 pcs.
SAAR646		1 g	-	-	667 PPM	53 PPM	-	100 pcs.
SAAR1647		1 g	-	-	1807 PPM	33 PPM	-	100 pcs.
SAAR1648		1 g	-	-	61 PPM	107 PPM	-	100 pcs.
SAAR1650		1 g	-	-	102 PPM	100 PPM	-	100 pcs.
SAAR1651	501-645	1 g	-	-	0.0210 %	0.0057 %	-	100 pcs.
SAAR1652		1 g	-	-	0.0530 %	0.0045 %	-	100 pcs.
SAAR1653		1 g	-	-	27 PPM	12 PPM	-	100 pcs.
SAAR1656		1 g	-	-	118 PPM	244 PPM	-	100 pcs.
SAAR654		0.5 g	-	-	534 PPM	65 PPM	-	100 pcs.
SAAR655	502-016 502-072	0.5 g	-	-	0.3930 %	0.0021 %	-	100 pcs.
SAAR656	501-552	0.5 g	-	-	0.0213 %	0.0081 %	-	100 pcs.
SAAR657		0.5 g	-	-	89 PPM	70 PPM	-	100 pcs.
SAAR659		0.5 g	-	-	66 PPM	120 PPM	-	100 pcs.
SAAR660	501-551	0.5 g	-	-	0.0060 %	0.0045 %	-	100 pcs.
SAAR661		0.5 g	-	-	167 PPM	9 PPM	-	100 pcs.
SAAR662		0.5 g	-	-	686 PPM	28 PPM	-	100 pcs.
SAAR663		0.5 g	-	-	0.1795 %	0.0042 %	-	100 pcs.
SAAR668		0.5 g	-	-	31 PPM	14 PPM	-	100 pcs.

Säntis No.	Leco No.	Size	C	H	N	O	S	Qty./ pack.
SAAR670		0.5 g	-	-	123 PPM	255 PPM	-	100 pcs.
SAAR675		0.5 g	-	-	0.2480 %	0.0770 %	-	100 pcs.
SAAR676		0.5 g	-	-	1214 PPM	329 PPM	-	100 pcs.

Zirconium Pins for Hydrogen & Nitrogen & Oxygen

SAAR640		0.1 g	-	14 PPM	50 PPM	1405 PPM	-	100 pcs.
---------	--	-------	---	--------	--------	----------	---	----------

Hydrogen in Steel Pins

SAAR546	762-747	1 g	-	1.1 PPM	-	-	-	100 pcs.
SAAR555	502-061	1 g	-	3.4 PPM	-	-	-	100 pcs.
SAAR556	502-060 501-529	1 g	-	5.0 PPM	-	-	-	100 pcs.
SAAR558		1 g	-	7.3 PPM	-	-	-	100 pcs.

Titanium Pins for Hydrogen, Oxygen, Nitrogen

SAAR637		0.1 g	-	12 PPM	101 PPM	509 PPM	-	100 pcs.
SAAR641		0.1 g	-	0.0039 %	0.0354 %	0.1910 %	-	100 pcs.
SAAR642		0.1 g	-	0.0050 %	0.0052 %	0.1140 %	-	100 pcs.
SAAR647		0.1 g	-	0.0090 %	0.0038 %	0.1750 %	-	100 pcs.
SAAR648	502-154	0.1 g	-	0.0161 %	0.0150 %	0.1580 %	-	100 pcs.
SAAR649		0.1 g	-	0.0200 %	0.0153 %	0.1610 %	-	100 pcs.
SAAR650	762-741	0.1 g	-	0.0061 %	0.0101 %	0.1555 %	-	100 pcs.
SAAR651	502-024	0.1 g	-	28 PPM	49 PPM	1320 PPM	-	100 pcs.

Oxygen in Copper Pins

SAAR147		1 g	-	-	-	10 PPM	11 PPM	100 pcs.
SAAR148	501-147	1 g	-	-	-	489 PPM	6 PPM	100 pcs.
SAAR149	501-148	1 g	-	-	-	323 PPM	11 PPM	100 pcs.

We can supply other standards with certificates of analysis upon request. So please don't hesitate to contact us.



COBRE

PARA ANALIZADORES ELEMENTALES

Listado de Artículos



EQUILAB

SANTIS
analytical

Article List for Copper for Elemental Analyzer

valid from 1st May 2010

Säntis No.	Description	Qty./ pack.
Copper Oxide		
SA99060100	Copper oxide wires 0,7 mm	Box 1 kg
SA99060101	Copper oxide wires 0,7 mm	Box 100 g
SA99060101B	Copper oxide wires 0,7 mm	Box 50 g
SA99060105	Copper oxide wires 0,7 mm	Box 250 g
SA99060300	Copper oxide wires 0,5 mm	Box 1 kg
SA99060301	Copper oxide wires 0,5 mm	Box 100 g
SA99060301B	Copper oxide wires 0,5 mm	Box 50 g
SA99060305	Copper oxide wires 0,5 mm	Box 250 g
SA99060306	Copper oxide granulat	Box 1 kg
SA99060307	Copper oxide granulat	Box 100 g
SA501170	Rare earth copper oxide	Vial 50 g

Reduced Copper (packed in Argon)			
SA02400117	Reduced copper granulate	Vial	681g
SA33835311	Electrolytic copper	Vial	90 g
SA33835311B	Electrolytic copper	Vial	70 g
SA33835312	High quality copper	Ampoule	50 g
SA33835305	High quality copper	Ampoule	220 g
SA33835335	High quality copper	Vial	250 g
SA99060102	Reduced copper wires 0,7 mm	Ampoule	100 g
SA99060103	Reduced copper wires 0,7 mm	Vial	250 g
SA99060104	Reduced copper wires 0,7 mm	Vial	200g
SA99060204	Reduced copper wires 0,7 mm (for sulfur analysis)	Ampoule	40 g
SA99060204B	Reduced copper wires 0,35 mm (for sulfur analysis)	Ampoule	40 g
SA99060205	Reduced copper wires 0,5 mm (for sulfur analysis)	Ampoule	40 g
SA99060302	Reduced copper wires 0,5 mm	Vial	100 g
SA99060302B	Reduced copper wires 0,5 mm	Ampoule	100 g
SA99060303	Reduced copper wires 0,5 mm	Vial	250 g
SA502304500E	Reduced Copper wires 0,5 mm	Vial	500 g
SA7710	Reduced Copper wires 0,5 mm	Vial	500 g
SA99060304	Reduced copper granulate	Vial	100g
SA99060304B	Reduced copper granulate	Ampoule	100g

Copper Turnings			
SA501621W	Copper Turnings Degassed	Box	200 g


We can manufacture other copper than those listed above. So please don't hesitate to contact us.



Perkin Elmer

Productos para ICP

Listado de Artículos



Säntis No.	Perkin Elmer No.	Description
SAPG30005	0047-2048	TORCH, ONE PIECE, with BALL JOINT
SAPG30012	0047-2047	SPRAY CHAMBER, GLASS SCOTT
SAPG30014		SPRAY CHAMBER, GLASS SCOTT, for conc. neb.
SAPG30015		SPRAY CHAMBER, TOP PORT ANGLED
SAPG30016		SPRAY CHAMBER, GLASS SCOTT
SAPG30018DV	N812-2188	CYCLONIC SPRAY CHAMBER
SAPG300-18DVJ		CYCLONIC SPRAY CHAMBER, 11.3 mm Exit, Water-Cooled
SAPG300-18DVQ		CYCLONIC SPRAY CHAMBER, 11.3 mm Exit, Quartz
SAPG30019		CYCLONIC SPRAY CHAMBER
SAPG30019J		CYCLONIC SPRAY CHAMBER, WATER COOLED
SAPG30019MS	N812-0150	CYCLONIC SPRAY CHAMBER
SAPG30019MSJ		CYCLONIC SPRAY CHAMBER, WATER COOLED
SAPG30019Q		QTZ CYCLONIC SPRAY CHAMBER
SAPG30019MSQ		QTZ CYCLONIC SPRAY CHAMBER
SAPG30019QSMJ		CYCLONIC SPRAY CHAMBER, QUARTZ
SAPG30021	0047-3052	TORCH EXTENSION
SAPG30022		TORCH EXTENSION with OPTICAL PORT
SAPG30023	0047-3560	TORCH EXTENSION, SLOTTED
SAPG30025	0047-2033	TORCH BONNET
SAPG30028		NEB. ADAPTER for CYCLONIC SPRAY CHAMBER
SAPG30029	N077-6052	CYCLONIC SPRAY CHAMBER, no BAFFLE with neb. adapter
SAPG30030	N077-6053	CYCLONIC SPRAY CHAMBER, BAFFLED with neb. adapter
SAPG30030Q		CYCLONIC SPRAY CHAMBER, BAFFLED, QUARTZ
SAPG30031	N077-6090	CYCLONIC SPRAY CHAMBER, 4mm BAFFLE
SAPG30032	N077-1501	CYCLONIC SPRAY CHAMBER, BAFFLED, 4300V
SAPG31003	N058-0528	DEMOUNTABLE TORCH TUBE, TYPE II
SAPG31003W	N068-1690	DEMOUNTABLE TORCH TUB., type II wear metals
SAPG31004	N069-0568	TORCH BODY, OPTIMA
SAPG31004W		TORCH TUBE FOR OPTIMA WEAR METALS
SAPG31005	N069-5295	TORCH BONNET, FLAT EDGE, Optimas bef. 9/94
SAPG31007	N058-0233	DEMOUNTABLE TORCH TUBE
SAPG31009	N058-0376	DEMOUNTABLE TORCH TUBE, TYPE I slotted
SAPG31011	N069-5379	TORCH BODY, OPTIMA xl
SAPG31011W		OPTIMA XL TORCH TUBE AWAR METALS
SAPG31012		OPTIMA DV TORCH TUBE 2 TOP SLOTS
SAPG31012W		310-12 WEAR METALS
SAPG31013	N069-5456	Torch Bonnet, Flat Edge (Optima post 9/94)
SAPG31014	N069-1662	TORCH BODY, DV 1 SLOT
SAPG31014W		OPTIMA DV WEAR METALS TORCH
SAPG310/3SL		OPTIMA DV TORCH, 3 SLOTS
SAPG31015	N069-1664	TORCH BONNET
SAPG31016	N077-0338	TORCH BODY, OPTIMA 2000/4000
SAPG31016W		TORCH BODY, OPTIMA 2000/4000 Wear Metals

Säntis No.	Perkin Elmer No.	Description
SAPG31017	N077-1500	TORCH BODY, RADIAL, 4300V
SAPG31018	N077-0344	TORCH BODY, SHORT VERSION
SAPG31019	N077-0343	TORCH BODY, NO SLOT
SAPG31020	N077-5297	COPPER FOIL for TORCH BODY OPTIMA 2000/4000
SAPG31024	N058-2235	QUARTZ BONNET
SAPG31025	N077-5289	TORCH BONNET
SAPG30029Q		CYCLONIC SPRAY CHAMBER, NO BAFFLE, QUARTZ
SAPG31035	N058-2248	INJECTOR SUPPORT ADAPTER w/O'RINGS
SAPG31036	N058-0539	O'RING KIT for INJECTOR SUPPORT ADAPTER
SAPG31037		INJECTOR HOLDER DV
SAPG31038	N069-0769	O-RING KIT FOR OPTIMA XL/DV
SAPG31039	N069-0780	INJECTOR SUPPORT ADAPTER w/O'RINGS
SAPG31040	N077-5227	ALUMINIA INJECTOR, 0.8mm
SAPG31041	N077-5228	ALUMINIA INJECTOR, 1.2mm
SAPG31042	N077-5177	ALUMINIA INJECTOR, 2.0mm
SAPG31043	N077-0437	O'RING KIT FOR TORCH
SAPG31051	N078-0130	Quartz Torch 1 slot Optima 8x00 – IB , ICP-OES Optima 8000
SAPG32030		QUARTZ INJECTOR TIP, 1.2mm TYPE I
SAPG32031	N058-0449	QUARTZ INJECTOR TIP, 1.6mm TYPE I
SAPG32032		QUARTZ SAMPLE INJECTOR TUBE 2.0 TYPE I
SAPG32033	N068-1631	QUARTZ INJECTOR TIP, 1.2mm
SAPG32034	0047-3292	QUARTZ INJECTOR TIP, 1.6mm
SAPG32035	N812-0529	QUARTZ INJECTOR, 2.0mm
SAPG32040	N069-5441	QUARTZ SAMPLE INJECTOR, 2.0mm
SAPG32041	N069-5442	QUARTZ SAMPLE INJECTOR, 2.0mm
SAPG32043	N069-5448	QUARTZ SAMPLE INJECTOR, 1.6mm
SAPG32044	N069-5449	QUARTZ SAMPLE INJECTOR, 1.6mm
SAPG32046	N069-5460	QUARTZ SAMPLE INJECTOR, 1.2mm
SAPG32048	N069-5461	QUARTZ INJECTOR 0.8 mm XL/DV
SAPG32050	N069-5362	ALUMINA INJECTOR for OPTIMA XL/DV
SAPG32051	N812-6039	ALUMINA INJECTOR for ELAN
SAPG32052	N812-6040	ALUMINA INJECTOR for ELAN
SAPG32053	N812-6041	ALUMINA INJECTOR for ELAN
SAPG32054	N077-1531	Alumina injector for Perkin Elmer 4300V
SAPG32060	N068-1234	PURGE EXTENSION
SAPG32062	0047-3187	PURGE TUBE
SAPG32064	0047-3561	PURGE TUBE FOR SLOTTED EXTENSION
SAPG32065		Purge Viewing Window - XL
SAPG32066	N058-1497	PURGE EXTENSION WINDOW
SAPG32067	N069-1646	PURGE TUBE, DV
SAPG32068	N069-1678	PURGE TUBE WINDOW, AXIAL
SAPG32069	N069-1689	PURGE TUBE WINDOW, RADIAL
SAPG32071	B081-0377	RADIAL PURGE WINDOW, 2000

Säntis No.	Perkin Elmer No.	Description
SAPG32073	N069-0672	PURGE TUBE WINDOW, RADIAL, SHORT
SAPG32074	N077-5014	QUARTZ INJECTOR, 2.0mm
SAPG32075	N077-5224	QUARTZ INJECTOR, 3.0mm
SAPG32076	N077-5223	QUARTZ INJECTOR, 1.6mm
SAPG32077	N077-5226	QUARTZ INJECTOR, 1.2mm
SAPG32078	N0775225	QUARTZ INJECTOR, 0.8mm
SAPG32080	N077-0322	RADIAL PURGE WINDOW, 4000
SAPG32081	N077-0944	RADIAL PURGE TUBE, FROSTED
SAPG40005	N810-2001	SHORT FORM TORCH 135mm overall length
SAPG40007	N810-2006	LONG FORM TORCH 151mm overall length
SAPG40009		LONG FORM TORCH 150mm overall length
SAPG40010	N812-2006	TORCH BODY, ELAN
SAPG40012	N810-2031	SCOTT SPRAY CHAMBER
SAPG40013		SCOTT SPRAY CHAMBER, QUARTZ
SAPG40014		SCOTT SPRAY CHAMBER, Water Cooled
SAPG40016		SCOTT SPRAY CHAMBER, Water Cooled, Quartz
SAPG40020	WE02-5221	CYCLONIC SPRAY CHAMBER, 6100DRC, QUARTZ
SAPG40020J		CYCLONIC SPRAY CHAMBER, 6100DRC, Jack, QUARTZ
SAPG40025		TORCH BONNET
SAPG40030	WE02-3948	INJECTOR, 6100DRC, 2.0mm
SAPG40031	WE02-7030	QUARTZ INJECTOR, 6100DRC & DRC II 0,85 m
SAPG40032	WE02-7005	QUARTZ INJECTOR, 6100DRC & DRC II 1,5 m
SAPG40035	N812-0116	INJECTOR SUPPORT/TORCH BASE
SAPG40036	N812-0100	SCIEX O-RING ADAPTER KIT
SAPG40037	WE023951	INJECTOR SUPPORT
SAPG40038		O'RING KIT FOR THE INJECTOR SUPPORT
SAPG40010HE	W100-8384	TORCH BODY, HIGH EFFICIENCY
SAPG50070P	0047-2022	C-Type CONCENTRIC NEBULIZER, 1 lpm Ar, 50psi, 1mL/min
SAPG50071	0047-2022	C-TYPE CONCENTRIC NEBULIZER, 1 lpm Ar, 30psi, 2mL/min
SAPG50072	0047-2022	CONCENTRIC NEB., 30psi, 1L/min Ar, 3mlpm
SAPG50080	N068-1574	CONCENTRIC NEB-, 30psi, 0.7 L/min Ar, 3mlpm

The use of Perkin Elmer™ part numbers is intended for convenience only, and does not imply that the products are Perkin Elmer™ origin. All products are guaranteed to be of high quality and suitable for use in the stated instrument. Perkin Elmer™ is a registered trade mark.